

# ADVIEW 2 SIMPLE INTEGRATION

## PURPOSE

The purpose of this document is to outline how to integrate your ADview 2 into our free, downloadable ADviewer software. In addition to this documentation we have a video walkthrough available the URL below:

<https://www.adctoday.com/adview2-simple-emr-integration>

## STEP 1: DOWNLOAD SOFTWARE


The first step is to download the ADviewer software from our website. You can download software via the link below:

<https://www.adctoday.com/adview2>

On this page head over to the "EMR" tab. Along the "Simple" EMR column scroll to the bottom and click "Download Now":


Recommended For	Attaching a PDF to the patient record.	Directly sending vitals to the specific fields in the patient record.	People who have HL7 translation software in use (ex Mirth)
Cost	Free	\$\$	\$\$\$\$
Complexity	Simple	Intermediate	Advanced
Pros	Quick and easy	True integration to your EMR fields.	Integrate into existing EMR IT infrastructure.
Cons	Simple file attachments leave data siloed in the patient record.	May involve costs from your EMR vendor	Requires a robust IT infrastructure which can be complex.
Support level you'll need in house	No support needed.	EMR Vendor support will be needed.	Onsite IT support and possible vendor support will be needed.
Timeline	1-2 hours max	1-2 weeks	Depends on IT infrastructure complexity and requirements.
Additional Equipment Needed	Wifi Dongle or Hardline Ethernet	Wifi Dongle or Hardline Ethernet and Barcode scanner (9005SCAN)	Wifi Dongle or Hardline Ethernet and Barcode scanner (9005SCAN)
	<a href="#">Download Now</a>	<a href="#">Notify Me</a>	Call (800)ADC-2270 for more info.

Once this is downloaded you can double right click to install. Once installed, the program will automatically run and you will see an icon on your desktop called "Adviewer":



## STEP 2: CONNECT THE ADVIEW 2 TO THE ADVIEWER SOFTWARE

After installation of the software it will automatically run on your computer. You will be presented with the screen below:



Right now, the key information we want to record is your systems IP address and port number which is displayed in the upper left corner of the screen:



This is the information we will need to load into the ADviews Advanced configuration screens.

With this information, we can now configure our Adviev to send readings to this address. The first step is to plug in the Adviev 2 to a PC configured to access the Advanced Configuration. Plug a USB cable to the Adviev 2's USB service port. This is marked with a 'Gear' icon.

If you do not have a PC configured to access the Advanced Configuration please refer to the link below for how to set this up:

<insert link>

Once connected we can access the Advanced Configuration by typing the following address in you web browsers address bar:

<http://model260.local>

You will be presented with the following logon:

## Login

The login form is titled "Enter username and password to proceed". It contains a "User" dropdown menu with "Medical" selected, a "Password" text input field with "Password" entered, and a "Sign in" button.

Drop down the user to the "Service" user and enter the password for the service account. The default password is 'service'.

Once logged in navigate to 'Device Configuration → EMR Settings':

The screenshot shows the SunTech Model 260 web interface. The top navigation bar includes "SunTech Model 260", "Overview", "Device Configuration", "Administration", and "Help", along with a "Logout" button. The "Device Configuration" menu is open, showing options like "Ethernet", "Wi-Fi", "NIBP calibration", "Measurements / Display", "EMR settings", "EMR server certificates", "EMR client certificates", "EMR connection test", "Date / Time and Language", and "Power management". A red arrow points to the "EMR settings" option. The main content area shows "Device information" with a "General" tab selected, displaying fields for "Serial number", "Ethernet address", "Uptime", "Battery charge", and "Kernel". A battery status bar at the bottom indicates "100% - Cycles: 5".


This page will allow us to enter the IP address and port number we recorded on our Adviewer screen. There are two sections we will want to load this information. First once under the 'Configure DEC' heading and another under the 'Configure PDQ' heading.:

The screenshot displays the configuration interface for SunTech Model 260. The top navigation bar includes 'SunTech Model 260', 'Overview', 'Device Configuration', 'Administration', and 'Help', along with a 'Logout' button. The main content area is divided into two sections:

- Configure DEC (Reading Reports):** This section contains several input fields:
  - IP/Hostname:** 192.168.200.155
  - Port:** 9999 (indicated by a red arrow)
  - Use SSL:**
  - Timeout:** 1
  - Retries:** 30
  - Retry interval:** 1
  - MSH-5 Application name:** ADCADVIEWER
  - MSH-6 Facility name:** ADCTRIAGE
  - Prefix for OBR-3.1:** Prefix for Filler Entity ID
  - OBR-3.2:** Filler Namespace ID
  - OBR-3.3:** Filler Universal ID
  - OBR-3.4:** Filler Universal ID type
- Configure PDQ (Patient Queries):** This section contains two input fields:
  - IP/Hostname:** 192.168.200.155
  - Port:** 9999 (indicated by a red arrow)

Once you have entered the IP and port information into both fields you can click 'Apply' in the bottom left of the screen.

Lastly, if you are not using a barcode scanner to scan patient IDs then you will want to navigate to 'Device Configuration → Measurements / Display' and unclick 'Use HL7 connect to validate scanned patient IDs':



We are all set. You can now log out of the advanced configuration and disconnect your USB cable.

### STEP 3: SENDING READINGS TO THE ADVIEWER SOFTWARE

With the Adviev 2 configured we can now send readings to the Adviev software. First thing we want to make sure of is that the Adviev is connected to the network either by a hardwired Ethernet cable OR using our optional WiFi dongle.

You can now go ahead and take a reading on the device. Once the reading is taken you can send the reading over the network by pressing the 'Memory' button then the 'center nob' on the front of the device.



You will see a flashing "EMR" icon with a check next to the it:



If successful you will see a check box and a small envelope on the screen:



If not successful you will see an error message displayed on the screen. Most errors are related to the Adviver being disconnected from the network, the advanced configuration is not pointing to the correct address, OR the Adviver software is not running.

At this point you will now see the reading on screen within the Adviver software.

You now have the ability to export the reading either to a local printer OR a PDF file which you can use to attach to the patients record.

That is it!