Diagnostix™ ADSTATION™

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5660T

Mode d'emploi

 ADC
55 Commerce Drive
Hauppauge, NY 11788

 mdi Europa GmbH
Langenhagener Str. 71
30855 Langenhagen
GERMANY

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tel: 631-273-9600
toll free: 1-800-232-2670
fax: 631-273-9659

www.adctoday.com
email: info@adctoday.com

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Telephone: 631-273-9600 • Fax: 631-273-9659

Email: info@adctoday.com

Recommended separation distances between portable and mobile RF communications equipment and the Diagnostic Station.			
The Diagnostic Station is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Diagnostic Station can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Diagnostic Station as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter	Separation distance according to frequency of transmitter (m)		
(W)	150 kHz to 80 MHz $a=12\sqrt{P}$	80 kHz to 800 MHz $a=1,2\sqrt{P}$	800 kHz to 2.7 GHz $a=2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations, Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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Guidance and manufacturer's declaration - Electromagnetic immunity			
The Diagnostic Station is intended for use in the electromagnetic environment specified below. The customer or the user of the Diagnostic station should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	Pass	Pass	<p>Portable and mobile RF communications equipment should not be used closer to any part of the Diagnostic Station, including the cables, than the recommended distance, which is calculated using the equation applicable to the transmitter frequency. Recommended separation distance</p> $d = 1.2 \sqrt{P} \text{ 150 KHz to 80 MHz}$ $d = 1.2 \sqrt{P} \text{ 80 KHz to 800 MHz}$ $d = 1.2 \sqrt{P} \text{ 800 KHz to 2,7 GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b interference may occur in the vicinity of equipment marked with the following symbol:</p> <p>Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and the recommended distance is given in meters (m). Field strengths from fixed RF transmitters determined by an electromagnetic site survey a, should be less than the compliance level in each frequency range b. Interference may occur in the vicinity of devices marked with the following symbol.</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m	
Proximity fields from RF wireless communications equipment	Pass	Pass	
<p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2: These guidelines may not apply in all situations.</p> <p>Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> <p>a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and landmobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Diagnostic Station is used exceeds the applicable RF compliance level above, the Diagnostic Station should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Diagnostic Station.</p> <p>b.) Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Guidance and manufacture's declaration - Electromagnetic emission		
The Diagnostic Station is intended for use in the electromagnetic environment specified below. The customer or the user of the Diagnostic station should assure that it is used in such an environment.		
Emission Test	Compliance	Electromagnetic environment - guidance
RF-Emission CISPR 11	Group 1	The Diagnostic Station uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF-Emission CISPR 11	Class B	The Diagnostic Station is intended for use in all establishments, including residential areas and those directly connected to a public supply network that also supplies buildings used for residential purposes.
Harmonic emissions IEC 61000-3-2	Pass	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Pass	

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Guidance and manufacture's declaration - Electromagnetic immunity			
The Diagnostic Station is intended for use in the electromagnetic environment specified below. The customer or the user of the Diagnostic station should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environmental guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Con: ± 8 kV Air: ± 15 kV	Con: ± 8 kV Air: ± 15 kV	Floors should be wood, concrete or ceramic tile. If floor is covered with synthetic material, the relative humidity should be at least 30%. The quality of the supply voltage should be that of a typical business or hospital environment.
Electrical fast transient/burst IEC 61000-4-4	5/50 ns, 100 kHz; ± 2 kV	5/50 ns, 100 kHz; ± 2 kV	
Surge IEC 61000-4-5	1.2/50 (8/20) μ s LIL: ± 1.0 kV LIG: ± 2.0 kV	1.2/50 (8/20) μ s LIL: ± 1.0 kV LIG: ± 2.0 kV	The quality of the supply voltage should be that of a typical business or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT for 0.5 cycle (1 phase) 0% UT for 1 cycle 70% UT for 25/30 cycles (50/60 Hz)	0% UT for 0.5 cycle (1 phase) 0% UT for 1 cycle 70% UT for 25/30 cycles (50/60 Hz)	The quality of the supply voltage should be that of a typical business or hospital environment.
	0% UT for 250/300 cycles (50/60 Hz)	0% UT for 250/300 cycles (50/60 Hz)	
Power frequency (50Hz/60Hz) magnetic field IEC 6000-4-8	30 A/m 50 Hz	30 A/m 50 Hz	Mains frequency magnetic fields should be at a level characteristic of a typical location in a typical commercial hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level. v			

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www.adctoday.com/support/warranty-registration.

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