

# DIAGNOSTICS

2006 FEATURE CHART



## Diagnostic Products Replacement Parts Matrix



<b>ADC Item#</b>	<b>5110</b> Pocket Otoscope/ Ophthalmoscope Set	<b>5111</b> Pocket Otoscope Set	<b>5112</b> Pocket Ophthalmoscope Set	<b>5210</b> Standard Diagnostic Set	<b>5211</b> Otoscope	<b>5215</b> Complete Diagnostic Set	<b>5312</b> Dermascope
<b>Battery Handle</b>	5160	5160	5160	5260	5260	5260	5260
<b>Ophthalmoscope Head</b>	5140	N/A	5140	5240	N/A	5240	N/A
<b>Otoscope Head</b>	5120	5120	N/A	5220	N/A	5220	N/A
<b>Dermascope Head</b>	N/A	N/A	N/A	N/A	N/A	N/A	5312-3
<b>WA Adaptor for Disp. Specula</b>	N/A	N/A	N/A	5211-AD	5211-AD	5211-AD	N/A
<b>Ear Specula, 2.5mm</b>	N/A	N/A	N/A	5211-1	5211-1	5211-1	N/A
<b>Ear Specula, 3.5mm</b>	N/A	N/A	N/A	5211-2	5211-2	5211-2	N/A
<b>Ear Specula, 4.5mm</b>	N/A	N/A	N/A	5211-3	5211-3	5211-3	N/A
<b>Bent Arm Illuminator</b>	N/A	N/A	N/A	N/A	N/A	5215-TH	N/A
<b>Metal Tongue Depressor</b>	N/A	N/A	N/A	N/A	N/A	5215-TD	N/A
<b>Laryngeal Mirror, 3mm</b>	N/A	N/A	N/A	N/A	N/A	5215-M3	N/A
<b>Laryngeal Mirror, 4mm</b>	N/A	N/A	N/A	N/A	N/A	5215-M4	N/A
<b>Nasal Specula</b>	N/A	N/A	N/A	N/A	N/A	5215-NS	N/A
<b>Tongue Depressor Attachment</b>	N/A	N/A	N/A	N/A	N/A	5215-TO	N/A
<b>Lamp, 2.5v for Ophthalmoscope</b>	4531-1	N/A	4531-1	5210-4	N/A	5210-4	N/A
<b>Lamp, 2.5v for Otoscope</b>	4530-1	4530-1	N/A	5211-4	5211-4	5211-4	N/A
<b>Lamp, 2.5v for Dermascope</b>	N/A	N/A	N/A	N/A	N/A	N/A	5312-4
<b>Dermascope Plain Contact Plate</b>	N/A	N/A	N/A	N/A	N/A	N/A	5312-1
<b>Dermascope Grad. Contact Plate</b>	N/A	N/A	N/A	N/A	N/A	N/A	5312-2
<b>Zipper Case</b>	N/A	N/A	N/A	5210-ZC	N/A	5215-ZC	5312-ZC
<b>Disp. Specula, 2.5mm, 1000/pkg</b>	5181	5181	N/A	5181	5181	5181	N/A
<b>Disp. Specula, 4.0mm, 1000/pkg</b>	5184	5184	N/A	5184	5184	5184	N/A
<b>Insufflator Bulb, Otoscope</b>	5122	5122	N/A	5122	5122	5122	N/A
<b>Insufflator Fitting, Otoscope</b>	5121	5121	N/A	N/A	5121	N/A	N/A

## Diagnostic Products Feature Chart



<b>ADC Item#</b>	<b>5110</b> Pocket Otoscope/ Ophthalmoscope Set	<b>5111</b> Pocket Otoscope Set	<b>5112</b> Pocket Ophthalmoscope Set	<b>5210</b> Standard Diagnostic Set	<b>5211</b> Otoscope	<b>5215</b> Complete Diagnostic Set	<b>5312</b> Dermascope
<b>Voltage</b>	2.5v	2.5v	2.5v	2.5v	2.5v	2.5v	2.5v
<b>Warranty</b>	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	Lifetime
<b>Head/Handle Mount</b>	Threaded Screw Type	Threaded Screw Type	Threaded Screw Type	Bayonet Locking	Bayonet Locking	Bayonet Locking	Bayonet Locking
<b>Battery Type</b>	2 "AA" (included)	2 "AA" (included)	2 "AA" (included)	2 "C"	2 "C"	2 "C"	2 "C"
<b>Power Switch</b>	Slide Switch	Slide Switch	Slide Switch	Rheostatic Switch	Rheostatic Switch	Rheostatic Switch	Rheostatic Switch
<b>Handle</b>	Plastic	Plastic	Plastic	Chrome-plated brass	Chrome-plated brass	Chrome-plated brass	Chrome-plated brass
<b>Insufflator Fitting</b>	Included	Included	N/A	Attached	N/A	Attached	N/A
<b>Illumination</b>	Fiber Optic	Fiber Optic	Fiber Optic	Incandescent	Incandescent	Incandescent	Fiber Optic
<b>Otoscope Magnification</b>	2.5x	2.5x	N/A	4x	4x	4x	N/A
<b>Ophthalmoscope Apertures</b>	5	N/A	5	1	1	1	N/A
<b>Ophthalmoscope Lenses</b>	17	N/A	17	25	25	25	N/A
<b>Instrument Heads</b>	2	1	1	2	1	3 with accessories	1
<b>Latex Content</b>	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free
<b>Case</b>	Hard Plastic	Hard Plastic	Hard Plastic	Semi-hard Zipper	Leather Zipper	Semi-hard Zipper	Nylon Zipper
<b>Case Pack</b>	20	N/A	N/A	3	N/A	3	6
<b>Catalog Page</b>	87	87	87	89	89	88	90

## Diagnostic Products Cross Reference Chart



ADC Item#	5110 Pocket Otoscope/ Ophthalmoscope Set	5111 Pocket Otoscope Set	5112 Pocket Ophthalmoscope Set	5210 Standard Diagnostic Set	5211 Otoscope	5215 Complete Diagnostic Set	5312 Dermascope
<b>Mabis</b>	20-312-020	20-310-020	20-3111-020	20-664-000	20-810-000	20-205-060	N/A
<b>Mabis</b>	N/A	20-860-0XXX <sup>1</sup>	N/A	N/A	20-660-000	N/A	N/A
<b>Omron</b>	R-6500	K-248-XXX <sup>1</sup>	N/A	N/A	R-3185	R-3400	N/A
<b>Heine</b>	D-873.11.052	D-851.10.052	D-852.10.052	A-279.10.118	B-181.10.118	A-139.10.118	K-161.10.118
<b>Welch Allyn</b>	928210	22821	12821	98510	24212	93510	47351

Note 1: Color selection available, ADC black only



American Diagnostic Corporation

55 Commerce Drive, Hauppauge, New York 11788

Tel: 631-273-9600 • Fax: 631-273-9659

Toll Free: 1-800-ADC-2670

[www.adctoday.com](http://www.adctoday.com)

