

DIGITAL BLOOD PRESSURE

2006 FEATURE CHART



Digital Blood Pressure Summary Chart



ADC Item #	6012	6013	6014	6015	6016
Inflation	Semi-Automatic	Automatic	Advanced Automatic	Automatic	Advanced Automatic
Deflation	Auto Valve	Elec. Automatic	Elec. Automatic	Elec. Automatic	Elec. Automatic
Type	Upper Arm	Upper Arm	Upper Arm	Wrist	Wrist
Smart Logic	No	Yes	Yes	Yes	Yes
Pulse	Yes	Yes	Yes	Yes	Yes
Date/Time	No	Yes	Yes	No	Yes
Memory	Last Reading	14	2 X 30	90	2 X 30
Memory Zone	1	1	2	1	2
Average Mode	No	No	Yes (MAM)	No	Yes (MAM)
Sizes Available	1	3	1	1	1
Limb Range	9"-13"	9"-13" (Adult)	9"-13"	5"-7 5/8"	5"-7 5/8"
Quick Start Button	No	Yes	No	Yes	Yes
Cuff Storage Compartment	No	Yes	Yes	N/A	N/A
Carry Case	No	No	Yes	Yes	Yes
Battery Type	4 "AA"	4 "AA"	4 "AA"	2 "AAA"	2 "AAA"
A/C Adapter Available	No	Yes	Yes	No	No
Printer Available	No	No	Yes	No	Yes
Latex Content	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free
Warranty	5 Year	5 Year	5 Year	5 Year	5 Year
Case Pack	6	6	6	24	12
Carton Pack	3	3	3	N/A	3
Catalog Page	28	28	29	30	30
BHS Validation	Yes	Yes	Yes	No	Yes
Replaces Old Unit	6006N	6008N	6009	6011	6011

Digital Blood Pressure Cross Reference Chart



ADC	6012	6013	6014 ^A	6015	6016 ^A
A & D	UA702	UA787	N/A	UB401	N/A
A & D	N/A	UA779	N/A	UB328	N/A
A & D	N/A	UA767	N/A	UB322	N/A
A & D	N/A	UA767 Plus	N/A	N/A	N/A
Lumiscop	1060	1080	N/A	1091	N/A
Lumiscop	1093	1085M	N/A	1092	N/A
Lumiscop	N/A	1094	N/A	N/A	N/A
Lumiscop	N/A	1095	N/A	N/A	N/A
Lumiscop	N/A	1097	N/A	N/A	N/A
Mabis	04-203-001	04-206-00	N/A	04-235-001	04-248-001
Mabis	04-263-001	04-244-00	N/A	61-238-001	04-255-001
Mabis	61-268-001	04-330-001 ²	N/A	N/A	N/A
Mabis	N/A	04-320-001 ²	N/A	N/A	N/A
Mabis	N/A	04-310-001	N/A	N/A	N/A
Mabis	N/A	61-315-001 ²	N/A	N/A	N/A
Mabis	N/A	61-305-001	N/A	N/A	N/A
Omron	HEM-432C	HEM-773AC ¹	HEM-705CP ³	HEM-609	HEM-637 ⁴
Omron	HEM-412C	HEM-711AC ¹	N/A	N/A	HEM-630
Omron	HEM-422CR	HEM-712C	N/A	HEM-629	N/A
Omron	N/A	HEM-757 ²	N/A	N/A	N/A
Prestige	HM-10	HM-20	N/A	HM-40	N/A
Prestige	N/A	HM-25	N/A	N/A	N/A
Samsung	SBM-600M	SBM-300	N/A	SBM-200 ⁶	N/A
Samsung	N/A	SBM-400F	N/A	SBM-200TA ⁷	N/A
Samsung	N/A	SBM-500S ⁵	N/A	N/A	N/A
Samsung	N/A	SBM-600F	N/A	N/A	N/A

- *Note 1 Includes AC Adapter, ADC adapter sold separately
- *Note 2 High speed model. ADC's all work at "regular" speed
- *Note 3 Includes Printer, ADC printer sold separately
- *Note 4 Positioning sensor. This feature not available from ADC
- *Note 5 Manually selectable inflation level (no fuzzy logic)
- *Note 6 Date and time stamp (ADC has more memories, no date/time stamp)
- *Note 7 Touch Screen Display

Note A Only ADC offers Average mode (averages 3 readings)

Digital Blood Pressure Model #6012 ADvantage™ Semi-Automatic Blood Pressure Monitor

Replaces 6006N

Features:	Benefits:
Automatic Off	Conserves battery life
Extra Large Display	For easy reading
Single Large Button	Simplifies use
Latex-Free	Reduces risk of allergic reactions
Functional Check	Checks all display functions & calibration self-test
Short Tone Indicates Starting	Patient is aware of measurement
Long Tone Indicates Finished	Patient is aware of completion
Auto-Valve	Simplifies measurement
Digital Inflation Instructions	Guides patient with cuff pressure
Digital Error Messages	Guides patient when there's a problem
Weighs Less Than 2lbs	Lightweight for travel
D-Ring Contoured Cuff	Easy, self-application
Last Reading Memory	Simplifies record keeping
8 Bit Microprocessor	Improved logic and accuracy
5 Year Warranty	Longer consumer protection
Includes Record Book	For convenient record keeping
BHS Validated	*Customer's Assurance of Accuracy*

BHS (British Hypertension Society) validation was issued to only a handful of manufacturers including ours.

It is equivalent to an endorsement by the AMA.



ADVANTAGE™ 6012
SEMI-AUTOMATIC BP MONITOR

Digital Blood Pressure Model #6013 ADvantage™ Automatic Blood Pressure Monitor

Replaces 6008N

Features:	Benefits:
Automatic Off	Conserves battery life
Smart Logic	No need to set systolic inflation level
Extra Large Display	For easy reading
Memory Stores Last 14 Readings	Simplifies record keeping
Date/Time Stamp	Keeps more relevant records
Latex-Free	Reduces risk of allergic reactions
Functional Check	Displays battery warning
5 Beeps Indicate Starting	Patient is aware of measurement
Long Tone Indicates Finished	Patient is aware of completion
Digital Error Messages	Guides patient when there's a problem
Available in Small and Large Adult	For small & large adult patients
Cuff Storage Compartment	Keeps unit compact for neatness
A/C Adapter Option	Conserves battery life
D-Ring Contoured Cuff	Easy, self-application
8 Bit Microprocessor	Improved logic and accuracy
5 Year Warranty	Longer consumer protection
Includes Record Book	For convenient record keeping
BHS Validated	*Customer's Assurance of Accuracy*

BHS (British Hypertension Society) validation was issued to only a handful of manufacturers including ours.

It is equivalent to an endorsement by the AMA.



ADVANTAGE™ 6013
AUTOMATIC BP MONITOR

Digital Blood Pressure Model #6014 ADvantage™ Advanced Automatic Blood Pressure Monitor

Replaces 6009

Features:

Benefits:

<i>Automatic Off</i>	Conserves battery life
<i>Smart Logic</i>	No need to set systolic inflation level
<i>Two Zone Memory</i>	Ideal for use with two patients
<i>Date/Time Stamp</i>	Keeps more relevant records
<i>Latex-Free</i>	Reduces risk of allergic reactions
<i>Functional Check</i>	Displays battery warning
<i>5 Beeps Indicate Starting</i>	Patient is aware of measurement
<i>Long Tone Indicates Finished</i>	Patient is aware of completion
<i>Digital Error Messages</i>	Guides patient when there's a problem
<i>Cuff Storage Compartment</i>	Keeps unit compact for neatness
<i>D-Ring Contoured Cuff</i>	Easy, self-application
<i>A/C Adapter Option</i>	Conserves battery life
<i>Extra Large Display</i>	For easy reading
<i>2 Zone 30 Reading Memory</i>	Simplifies record keeping
<i>Printer Option</i>	For easy record keeping
<i>Average Mode</i>	Improved accuracy
<i>Large Adult Replacement Cuff Available</i>	For large adult patients
<i>8 Bit Microprocessor</i>	Improved logic and accuracy
<i>5 Year Warranty</i>	Longer consumer protection
<i>Includes Record Book</i>	For convenient record keeping
<i>*BHS Validated*</i>	<i>*Customer's Assurance of Accuracy*</i>

BHS (British Hypertension Society) validation was issued to only a handful of manufacturers including ours.

It is equivalent to an endorsement by the AMA.



ADVANTAGE™ 6014P
Optional thermal Printer



ADVANTAGE™ 6014
AUTOMATIC BP MONITOR

Digital Blood Pressure Model #6015
ADvantage™ Basic Wrist Blood Pressure Monitor

Replaces 6011

Features:

Automatic Off

Smart Logic

Memory Stores Last 90 Readings

Latex-Free

Functional Check

Digital Error Messages

Weighs Less Than 1lb

Includes Carry Case

5 Year Warranty

Includes Record Book

Single Button Operation

Measures at Wrist

Contoured Cuff

Benefits:

Conserves battery life

No need to set systolic inflation level

Simplifies record keeping

Reduces risk of allergic reactions

Displays battery warning

Guides patient when there's a problem

Compact & lightweight for travel

Convenient for travel

Longer consumer protection

For convenient record keeping

Simplifies use

Improved portability & convenience

Easy, self-application



ADVANTAGE™ 6015
AUTOMATIC BP MONITOR

**Digital Blood Pressure Model #6016
ADvantage™ Advanced Wrist Blood Pressure Monitor**

Replaces 6011

Features:

Benefits:

<i>Automatic Off</i>	Conserves battery life
<i>Smart Logic</i>	No need to set systolic inflation level
<i>Two Zone Memory</i>	Ideal for use with two patients
<i>Date/Time Stamp</i>	Keeps more relevant records
<i>Latex-Free</i>	Reduces risk of allergic reactions
<i>Long Tone Indicates Finished</i>	Patient is aware of completion
<i>Digital Error Messages</i>	Guides patient when there's a problem
<i>Contoured Cuff</i>	Easy, self-application
<i>Memory Stores Last 30 Readings Each</i>	Simplifies record keeping
<i>Printer Option</i>	For easy record keeping
<i>Average Mode</i>	Improved accuracy
<i>8 Bit Microprocessor</i>	Improved logic and accuracy
<i>5 Year Warranty</i>	Longer consumer protection
<i>Includes Record Book</i>	For convenient record keeping
<i>Weighs Less Than 1lb</i>	Compact & lightweight for travel
<i>Includes Carry Case</i>	Convenient for travel
<i>Measures at Wrist</i>	Improved portability & convenience
<i>*BHS Validated*</i>	<i>*Customer's Assurance of Accuracy*</i>

BHS (British Hypertension Society) validation was issued to only a handful of manufacturers including ours.

It is equivalent to an endorsement by the AMA.



**ADvantage™ 6016
AUTOMATIC BP MONITOR**



American Diagnostic Corporation

55 Commerce Drive, Hauppauge, New York 11788
tel: 631-273-9600 • fax: 631-273-9659
toll free: 1-800-ADC-2670

www.adctoday.com

