

Digital Blood Pressure Model #6015  
ADvantage™ Basic Wrist Blood Pressure Monitor



**Features:**

- Automatic Off
- Smart Logic
- Memory Stores Last 30 Readings with Date/Time Stamp
- Latex-Free
- Functional Check
- Digital Error Messages
- Weighs Less Than 1lb
- Includes Carry Case
- 5 Year Warranty
- Includes Record Book
- Single Button Operation
- Measures at Wrist
- Contoured Cuff

**Benefits:**

- Conserves battery life
- No need to set systolic inflation level
- Simplifies record keeping
- Reduces risk of allergic reactions
- Displays battery warning
- Guides patient when there's a problem
- Compact & lightweight for travel
- Convenient for travel
- Longer consumer protection
- For convenient record keeping
- Simplifies use
- Improved portability & convenience
- Easy, self-application

Digital Blood Pressure Model #6016  
ADvantage™ Advanced Wrist Blood Pressure Monitor



**Features:**

- Automatic Off
- Smart Logic
- Two Zone Memory
- Date/Time Stamp
- Latex-Free
- Long Tone Indicates Finished
- Digital Error Messages
- Contoured Cuff
- Memory Stores Last 30 Readings Each
- PC Software Available
- Average Mode
- 8 Bit Microprocessor
- 5 Year Warranty
- Includes Record Book
- Weighs Less Than 1lb
- Includes Carry Case
- Measures at Wrist
- \*BHS Validated\*

**Benefits:**

- Conserves battery life
- No need to set systolic inflation level
- Ideal for use with two patients
- Keeps more relevant records
- Reduces risk of allergic reactions
- Patient is aware of completion
- Guides patient when there's a problem
- Easy, self-application
- Simplifies record keeping
- For easy data analysis
- Improved accuracy
- Longer consumer protection
- For convenient record keeping
- Compact & lightweight for travel
- Convenient for travel
- Improved portability & convenience
- \*Customer's Assurance of Accuracy\*

BHS (British Hypertension Society) validation was issued to only a handful of manufacturers including ours. Equivalent endorsement to the AMA.



ADC  
55 Commerce Drive  
Hauppauge, NY 11788



ADC (UK) Ltd.  
Unit 6, PO14 1TH  
United Kingdom

tel: 631-273-9600  
toll free: 1-800-232-2670  
fax: 631-273-9659

www.adctoday.com  
info@adctoday.com



6 3 4 7 8 2 0 5 4 4 7 2  
digbpf2004 - rev 3

# DIGITAL BLOOD PRESSURE

## FEATURE CHART III



## Digital Blood Pressure Summary Chart



ADC Item #	6012	6015	6016	6021	6022	6023
<b>Inflation</b>	Semi-Auto	Automatic	Advanced Auto	Automatic	Advanced Auto	Advanced Auto
<b>Deflation</b>	Auto Valve	Elec. Auto	Elec. Auto	Elec. Auto	Elec. Auto	Elec. Auto
<b>Type</b>	Upper Arm	Wrist	Wrist	Upper Arm	Upper Arm	Upper Arm
<b>Smart Logic</b>	No	Yes	Yes	Yes	Yes	Yes
<b>Pulse</b>	Yes	Yes	Yes	Yes	Yes	Yes
<b>Date/Time</b>	No	Yes	Yes	Yes	Yes	Yes
<b>Memory</b>	Last Reading	30	2 X 30	99	2 X 99	2 X 99
<b>Memory Zone</b>	1	1	2	1	2	2
<b>Average Mode</b>	No	No	Yes (MAM)	No	Yes (MAM)	Yes (MAM)
<b>Sizes Available</b>	1	1	1	3	1	2
<b>Limb Range</b>	9"-13" Adult	5"-7 5/8" Wrist	5"-7 5/8" Wrist	7.5"-9" Sm. Adult 9"-13" Adult 13"-17" Lg. Adult	9"-13" Adult	9"-13" Adult 13"-17" Lg. Adult
<b>Quick Start Button</b>	Yes	Yes	Yes	Yes	Yes	Yes
<b>Cuff Storage Compartment</b>	No	N/A	N/A	Yes	Yes	No
<b>Carry Case</b>	No	Yes	Yes	Yes	Yes	Yes
<b>Battery Type</b>	4 "AA"	2 "AAA"	2 "AAA"	4 "AA"	4 "AA"	4 "AA"
<b>A/C Adapter Available</b>	No	No	No	Yes	Yes	Yes
<b>Latex Content</b>	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free
<b>Warranty</b>	5 Year	5 Year	5 Year	5 Year	5 Year	5 Year
<b>Case Pack</b>	6	24	12	6	6	6
<b>Carton Pack</b>	3	12	3	3	3	3
<b>Catalog Page</b>	36	36	36	38	38	37
<b>BHS Validation</b>	Yes	No	Yes	Yes	Yes	Yes
<b>Replaces Old Unit</b>	6006N	6011	6011	6013	6014	6017
<b>PC Connectivity</b>	N/A	N/A	Yes	N/A	Yes	Yes

## Digital Blood Pressure Model #6023 ADvantage™ **ULTRA** Advanced Blood Pressure Monitor

### Features:

**Automatic Off**

**Smart Logic**

**Traffic Light Indicator**

**Date/Time Stamp**

**Latex-Free**

**Irregular Heartbeat Detection**

**Digital Error Messages**

**Wide Range Rigid Cuff**

**Advanced 2 Zone Memory**

**MAM Advantage™ Average Mode**

**PC Link Software**

**5 Year Warranty**

**AC Adapter**

**USB Cable**

**Deluxe Carry Case**

**\*BHS Validated**

### Benefits:

Conserves battery life

Automatically determines ideal inflation levels

Compares reading with National Standards

Keeps more relevant records

Reduces risk of allergic reactions

Signals appearance of irregular heartbeat

Guides patient when there's a problem

Easy, self-application for medium and large adults

Recalls last 99 readings for two people

Improved accuracy

Connects monitor to your PC to download and print data, graphs, charts, and tables from readings

Longer consumer protection

Conserves battery life

For easy connection to your PC

Convenient for travel and storage

Customer's Assurance of Accuracy

**\*BHS (British Hypertension Society) validation was issued to only a handful of manufacturers including ours.**

*It is equivalent to an endorsement by the AMA.*

ADVANTAGE™ 6023 **ULTRA**  
AUTOMATIC BP MONITOR





Digital Blood Pressure Model #6012  
ADvantage™ Semi-Automatic Blood Pressure Monitor

Features:	Benefits:
<b>Automatic Off</b>	Conserves battery life
<b>Extra Large Display</b>	For easy reading
<b>Single Large Button</b>	Simplifies use
<b>Latex-Free</b>	Reduces risk of allergic reactions
<b>Functional Check</b>	Checks all display functions & calibration self-test
<b>Short Tone Indicates Starting</b>	Patient is aware of measurement
<b>Long Tone Indicates Finished</b>	Patient is aware of completion
<b>Auto-Valve</b>	Simplifies measurement
<b>Digital Inflation Instructions</b>	Guides patient with cuff pressure
<b>Digital Error Messages</b>	Guides patient when there's a problem
<b>Weighs Less Than 2lbs</b>	Lightweight for travel
<b>D-Ring Contoured Cuff</b>	Easy, self-application
<b>Last Reading Memory</b>	Simplifies record keeping
<b>8 Bit Microprocessor</b>	Improved logic and accuracy
<b>5 Year Warranty</b>	Longer consumer protection
<b>Includes Record Book</b>	For convenient record keeping
<b>*BHS Validated*</b>	*Customer's Assurance of Accuracy*

**BHS (British Hypertension Society)** validation was issued to only a handful of manufacturers including ours.

It is equivalent to an endorsement by the AMA.



ADVANTAGE™ 6012  
SEMI-AUTOMATIC BP MONITOR

Digital Blood Pressure Model #6021  
ADvantage™ Automatic Blood Pressure Monitor

Features:	Benefits:
<b>Automatic Off</b>	Conserves battery life
<b>Smart Logic</b>	No need to set systolic inflation level
<b>Extra Large Display</b>	For easy reading
<b>Memory Stores Last 99 Readings</b>	Simplifies record keeping
<b>Date/Time Stamp</b>	Keeps more relevant records
<b>Latex-Free</b>	Reduces risk of allergic reactions
<b>Irregular Heartbeat Detection</b>	Signals appearance of irregular heartbeat
<b>Functional Check</b>	Displays battery warning
<b>5 Beeps Indicate Starting</b>	Patient is aware of measurement
<b>Long Tone Indicates Finished</b>	Patient is aware of completion
<b>Digital Error Messages</b>	Guides patient when there's a problem
<b>Available in Small and Large Adult</b>	For small & large adult patients
<b>Cuff Storage Compartment</b>	Keeps unit compact for neatness
<b>A/C Adapter Option</b>	Conserves battery life
<b>D-Ring Contoured Cuff</b>	Easy, self-application
<b>8 Bit Microprocessor</b>	Improved logic and accuracy
<b>5 Year Warranty</b>	Longer consumer protection
<b>Includes Record Book</b>	For convenient record keeping
<b>Deluxe Carry Case</b>	For convenient travel and storage
<b>*BHS Validated*</b>	*Customer's Assurance of Accuracy*

**BHS (British Hypertension Society)** validation was issued to only a handful of manufacturers including ours.

It is equivalent to an endorsement by the AMA.



ADVANTAGE™ 6021  
AUTOMATIC BP MONITOR