



ADLITE™

Diagnostic Penlights



REUSABLE PENLIGHTS

Three distinctly different models to suit your personal preferences.



ADLITE PRO™

Our finest pocket light, the ADLITE™ Pro is the physician's choice for evaluation of the pupillary reflex.

- LED lamp for bright, white illumination for true tissue rendering
- LED lamp provides 10,000 hour use eliminating the high cost of replacement lamps
- Durable CNC machined brass housing in a striking black matte finish
- Positive feel power switch
- Measures 5 5/8" long, 1/2" diameter
- Convenient Pocket Clip
- Complete with 2 "AAA" batteries and an attractive gift box
- Lifetime warranty on lamp cartridge, 2 year warranty on penlight

ORDERING INFORMATION

#355BK	Adlite Pro™, Physician's Penlight
--------	-----------------------------------



METALITE™

- Brushed aluminum housing
- Locking push toggle activated
- Measures 5 1/4" long by 1/2" diameter
- Two "AAA" batteries included

ORDERING INFORMATION

#352	METALITE™ Reusable Penlight
#352CYL	METALITE™ Reusable Penlight, 22/cylinder



ADLITE II™

- Plastic housing for lightweight comfort
- Pocket clip activated
- Measures 5" long by 1/2" diameter
- Two "AAA" batteries included

ORDERING INFORMATION

#354	ADLITE II™ Reusable Penlight
#354CYL	ADLITE II™ Reusable Penlight, 18/cylinder



ADLITE™

- Sealed concentric beam, disposable diagnostic penlight
- Activated by depressing pocket clip
- Available with or without pupil gauge
- Measures 5" long, 1/2" diameter
- Sold in 6 pack



ORDERING INFORMATION

#351	ADLITE™ Disposable Penlight, 6/pk
#351P	ADLITE™ Disposable Penlight with Pupil Gauge, 6/pk

ADLITE™ PLUS

- Sealed concentric beam, disposable diagnostic penlight
- Activated by depressing pocket clip
- White available with or without pupil gauge
- Measures 5" long, 1/2" diameter
- Sold individually or in display cylinder of 22
- Available in 8 colors



ORDERING INFORMATION

#356*	ADLITE™ Disp. Penlight, specify color
Specify: Black - BK, White - W, White w/Pupil Gauge - WP, Orange - OR, Purple - V, Red - R, Teal - TL, Royal Blue -RB	
#358AST	Disp. Penlight, 22/Cylinder - Assorted
#358TVB	Disp. Penlight, 22/Cylinder - Black, Purple, Teal
#358W	Disp. Penlight, 22/Cylinder - White
#358WP	Disp. Penlight, 22/Cylinder - White W/Pupil Gauge



55 Commerce Drive
Hauppauge, New York 11788
tel: 1-800-232-2670
fax: 631-273-9659

www.adctoday.com



COLOR CHOICES

BK	WP	OR	V	R	TL	RB
----	----	----	---	---	----	----