


DIAGNOSTICS

FEATURE CHART III



Diagnostic Products Replacement Parts Matrix



ADC Item#	5110E Pocket Economy Diagnostic Set	5110N Pocket Otoscope /Ophthalmoscope Set	5111N Pocket Otoscope Set	5112N Pocket Ophthalmoscope Set	5210 Standard Diagnostic	5211 Otoscope Set	5215 Complete Diagnostic	5312 Dermascope
Battery Handle	5160N	5160N	5160N	5160N	5211-5	5211-5	5211-5	5260
Ophthalmoscope Head	5140N	5140N	N/A	5140N	5240	N/A	5240	N/A
Otoscope Head	5120N	5120N	5120N	N/A	5220	5220	5220	N/A
Dermascope Head	N/A	N/A	N/A	N/A	N/A	N/A	N/A	5312-3
WA Adaptor for Disp. Specula	N/A	N/A	N/A	N/A	5211-AD	5211-AD	5211-AD	N/A
Ear Specula, 2.5mm	N/A	N/A	N/A	N/A	5211-1	5211-1	5211-1	N/A
Ear Specula, 3.5mm	N/A	N/A	N/A	N/A	5211-2	5211-2	5211-2	N/A
Ear Specula, 4.5mm	N/A	N/A	N/A	N/A	5211-3	5211-3	5211-3	N/A
Bent Arm Illuminator	N/A	N/A	N/A	N/A	N/A	N/A	5215-TH	N/A
Metal Tongue Depressor	N/A	N/A	N/A	N/A	N/A	N/A	5215-TD	N/A
Laryngeal Mirror, 3mm	N/A	N/A	N/A	N/A	N/A	N/A	5215-M3	N/A
Laryngeal Mirror, 4mm	N/A	N/A	N/A	N/A	N/A	N/A	5215-M4	N/A
Nasal Specula	N/A	N/A	N/A	N/A	N/A	N/A	5215-NS	N/A
Tongue Depressor Attachment	N/A	N/A	N/A	N/A	N/A	N/A	5215-T0	N/A
Lamp, 2.5v for Ophthalmoscope	5112N-4	5112N-4	N/A	5112N-4	5210-4	N/A	5210-4	N/A
Lamp, 2.5v for Otoscope	5111N-4	5111N-4	5111N-4	N/A	5211-4	5211-4	5211-4	N/A
Lamp, 2.5v for Dermascope	N/A	N/A	N/A	N/A	N/A	N/A	N/A	5312-4
Dermascope Plain Contact Plate	N/A	N/A	N/A	N/A	N/A	N/A	N/A	5312-1
Dermascope Grad. Contact Plate	N/A	N/A	N/A	N/A	N/A	N/A	N/A	5312-2
Case	925110E-00	925110N-00	925111N-00	925112N-00	5210-ZC	N/A	5215-ZC	5312-ZC
Disp. Specula, 2.75mm, 850/pkg	5182	5182	5182	N/A	5182	5182	5182	N/A
Disp. Specula, 4.25mm, 850/pkg	5185	5185	5185	N/A	5185	5185	5185	N/A
Insufflator Bulb, Otoscope	5122N	5122N	5122N	N/A	5122N	5122N	5122N	N/A
Insufflator Fitting, Otoscope	5121N	5121N	5121N	N/A	N/A	5121N	N/A	N/A

Diagnostic Products Feature Chart



ADC Item#	5110E Economy Diagnostic Set	5110N Pocket Otoscope/ Ophthalmoscope Set	5111N Pocket Otoscope Set	5112N Pocket Ophthalmoscope Set	5210 Standard Diagnostic	5211 Otoscope	5215 Complete Diagnostic Set	5312 Dermoscope
Voltage	2.5v	2.5v	2.5v	2.5v	2.5v	2.5v	2.5v	2.5v
Warranty	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	Lifetime
Head/Handle Mount	Threaded Screw Type	Threaded Screw Type	Threaded Screw Type	Threaded Screw Type	Bayonet Locking	Bayonet Locking	Bayonet Locking	Bayonet Locking
Battery Type	2 "AA" (included)	4 "AA" (included)	2 "AA" (included)	2 "AA" (included)	2 "C"	2 "C"	2 "C"	2 "C"
Power Switch	Slide Switch	Slide Switch	Slide Switch	Slide Switch	Rheostatic Switch	Rheostatic Switch	Rheostatic Switch	Rheostatic Switch
Handle	Plastic	Plastic	Plastic	Plastic	Chrome-plated brass	Chrome-plated brass	Chrome-plated brass	Chrome-plated brass
Insufflator Fitting	Sold Separately	Sold Separately	Sold Separately	Sold Separately	Attached	N/A	Attached	N/A
Illumination	Fiber Optic	Fiber Optic	Fiber Optic	Fiber Optic	Incandescent	Incandescent	Incandescent	Fiber Optic
Otoscope Magnification	2.5x	2.5x	2.5x	N/A	4x	4x	4x	N/A
Ophthalmoscope Aperatures	5	5	N/A	5	1	1	1	N/A
Ophthalmoscope Lenses	19	19	N/A	19	24	24	24	N/A
Instrument Heads	2	2	1	1	2	1	3 with accessories	1
Latex Content	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free
Case	Hard Plastic	Hard Plastic	Hard Plastic	Hard Plastic	Semi-hard Zipper	Leather Zipper	Semi-hard Zipper	Nylon Zipper
Case Pack	12	12	12	12	3	N/A	3	6

Diagnostic Products Cross Reference Chart



ADC Item#	5110E Pocket Economy Diagnostic Set	5110N Pocket Otoscope/ Ophthalmoscope Set	5111N Pocket Otoscope Set	5112N Pocket Ophthalmoscope Set	5210 Standard Diagnostic Set	5211 Otoscope	5215 Complete Diagnostic Set	5312 Dermascope
Mabis	N/A	20-312-020	20-310-020	20-311-020	20-664-000	20-810-000	20-205-060	N/A
Mabis	N/A	20-332-020	20-860-0XXX¹	N/A	N/A	20-660-000	N/A	N/A
Omron	N/A	R-6500	K-248-XXX¹	N/A	N/A	R-3185	R-3400	N/A
Heine	N/A	D-873.21.021 D-294.10.051	D-851.20.021 D-195.10.051	D-852.10.021 D-001.71.120 D-001.83.120	A-279.10.118 A-177.10.118 A-181.10.118	B-141.10.118 B-181.10.118	A-139.10.118	K-161.10.118 D-001.78.107
Welch Allyn	N/A	92820 95001	22841 25020	04900-U 12851	98510	20250 24212	93510	47351

Note 1: Color selection available, ADC black only



ADC
55 Commerce Drive
Hauppauge, NY 11788

tel: 631-273-9600
toll free: 1-800-232-2670
fax: 631-273-9659

www.adctoday.com

info@adctoday.com



ADC (UK) Ltd.
Unit 6, PO14 1TH
United Kingdom

