


LARYNGOSCOPES



Fiberoptic



MACINTOSH BLADES

- Stainless steel construction with Satin non-glare finish
- Fiberoptic light transmission provides cool, obstruction free viewing
- 2.5v Halogen illumination for true tissue color & long lasting lamp performance
- Removable light pipe facilitates cleaning and service
- Color coded green to signify compatibility ONLY with fiberoptic handles



ORDERING INFORMATION

Item#	Blade	Description	Dim A	Dim B
4070F	0	Preemie	80	55
4071F	1	Infant	92	67
4072F	2	Child	108	83
4073F	3	Med. Adult	126	101
4074F	4	Lrg. Adult	155	134

(Dimensions in millimeters)


MILLER BLADES

- Stainless steel construction with Satin non-glare finish
- Fiberoptic light transmission provides cool, obstruction free viewing
- 2.5v Halogen illumination for true tissue color & long lasting lamp performance
- Removable light pipe facilitates cleaning and service
- Color coded green to signify compatibility ONLY with fiberoptic handles

ORDERING INFORMATION

Item#	Blade	Description	Dim A	Dim B
4080F	0	Preemie	78	53
4081F	1	Infant	102	77
4082F	2	Child	155	103
4083F	3	Med. Adult	195	170
4084F	4	Lrg. Adult	205	180

(Dimensions in millimeters)



BATTERY HANDLES

- Satin chrome-plated brass handle with a knurled finish for a positive grip
- Light source in handle simplifies blade cleaning, and provides a brighter, whiter light
- Color coded green to signify compatibility ONLY with fiberoptic blades
- Available in traditional "C", "AA", and stubby ("AA")
- Batteries not included




ORDERING INFORMATION

Item#	Description
4065F	Medium "C" size Fiberoptic battery handle
4066F	Small "AA" size Fiberoptic battery handle
4067F	Stubby "AA" size Fiberoptic battery handle
4532-1	Fiberoptic Replacement Lamp, Pkg. of 1



Standard



MACINTOSH STANDARD

- Stainless steel construction with Satin non-glare finish
- The most popular of the curved blade designs

ORDERING INFORMATION

Item#	Blade	Description	Dim A	Dim B	Lamp
4070	0	Preemie	80	58	#4500
4071	1	Infant	87	65	#4500
4072	2	Child	108	86	#4500
4073	3	Med. Adult	128	105	#4501
4074	4	Lrg. Adult	159	137	#4501

(Dimensions in millimeters)


MILLER STANDARD

- Stainless steel construction with Satin, non-glare finish
- The most popular of the straight blade designs

ORDERING INFORMATION

Item#	Blade	Description	Dim A	Dim B	Lamp
4080-0	00	Small Preemie	65	43	#4500
4080	0	Preemie	75	53	#4500
4081	1	Infant	102	80	#4500
4082	2	Child	155	133	#4501
4083	3	Med. Adult	195	173	#4501
4084	4	Lrg. Adult	205	183	#4501

(Dimensions in millimeters)




WISCONSIN STANDARD

- Stainless steel construction with Satin, non-glare finish

ORDERING INFORMATION

Item#	Blade	Description	Dim A	Dim B	Lamp
4090	0	Preemie	77	55	#4500
4091	1	Infant	102	80	#4500
40915	1.5	Small Child	115	93	#4500
4092	2	Child	135	112	#4501
4093	3	Med. Adult	162	139	#4501
4094	4	Lrg. Adult	199	177	#4501

(Dimensions in millimeters)



BATTERY HANDLES

- Satin chrome-plated brass with a knurled finish for a positive grip
- Available in traditional "C", "AA", and stubby ("AA")
- Batteries not included



ORDERING INFORMATION

Item#	Description
4065	Medium "C" size Battery handle
4066	Small "AA" size Battery handle
4067	Stubby "AA" size Battery handle



ALL ADC® LARYNGOSCOPE BLADES & HANDLES CONFORM TO ASTM F965 AND ISO 7376-1 STANDARDS

Disposable



MACINTOSH

- Polycarbonate construction
- The most popular of the curved blade designs
- Works with our fiberoptic handles and any other conforming to ASTM-F965 or ISO 7376-1 standards
- Individually packaged in sealed pouch
- 20/carton, 120/case
- Non-sterile



ORDERING INFORMATION					
Item#	Blade	Description	Dim A	Dim B	Pkg
4072D-20	2	Child	114.3	89	20
4073D-20	3	Med. Adult	127	105	20
4074D-20	4	Lrg. Adult	155.5	133	20


(Dimensions in millimeters)

MILLER

- Polycarbonate construction
- The most popular of the straight blade designs
- Works with our fiberoptic handles and any other conforming to ASTM-F965 or ISO 7376-1 standards
- Individually packaged in sealed pouch
- 20/carton, 120/case
- Non-sterile

ORDERING INFORMATION					
Item#	Blade	Description	Dim A	Dim B	Pkg
4080D-20	0	Preemie	79.3	54	20
4081D-20	1	Infant	104.7	80	20
4082D-20	2	Child	152.4	133	20
4083D-20	3	Adult	196.8	171	20

(Dimensions in millimeters)



LED DISPOSABLE BATTERY HANDLE

- LED illuminated for bright, cool light
- Ribbed surface for positive grip
- Sized like a traditional "C" cell handle
- Powered by 3 "AAA" batteries (included)
- May be ETO sterilized
- Each in hermetically sealed pouch, 20/carton
- The lowest cost handle on the market



ORDERING INFORMATION	
Item#	Description
4065D-20	LED Disposable Battery Handle



FIBEROPTIC and STANDARD LARYNGOSCOPE SETS

Choose from either a 4 blade Macintosh, or 5 blade Miller set with both "AA" and "C" battery handles.

The Fiberoptic Laryngoscope Set features:

- 4 Macintosh or 5 Miller blades
- Both "AA" and "C" cell battery handles
- Foam lined fitted case with carry handle
- Batteries included

The Standard Laryngoscope Set features:

- 4 Macintosh or 5 Miller blades
- Both "AA" and "C" cell battery handles
- Foam lined fitted case with carry handle
- Batteries included



MACINTOSH and MILLER DISPOSABLE LARYNGOSCOPE SETS

Choose from either a 6 blade Macintosh or 8 blade Miller set with "C" battery handles.

The 4089D Fiberoptic Miller Disposable

Laryngoscope Set features:

- 8 Disposable Miller blades - 2 of each size (Sizes 0-3)
- Disposable battery handle, includes (3) "AAA" batteries
- Zippered carry case holds both blades and handle

The 4079D Fiberoptic Macintosh Disposable

Laryngoscope Set features:

- 6 Disposable Macintosh blades - 2 of each size (Sizes 2-4)
- Disposable battery handle, includes (3) "AAA" batteries
- Zippered carry case holds both blades and handle



ORDERING INFORMATION

#4079	Laryngoscope Set, Macintosh, includes one each of Mac #1, #2, #3, #4 blades; "AA" handle (#4066) & "C" handle (#4065) in hard case
#4089	Laryngoscope Set, Miller, includes one each of Miller #0, #1, #2, #3, #4, blades; "AA" handle (#4066), & "C" handle (#4065) in hard case
#4079F	Fiberoptic Laryngoscope Set, Macintosh, includes one each of Mac #1, #2, #3, #4 blades; "AA" handle (#4066F) & "C" handle (#4065F) in hard case
#4089F	Fiberoptic Laryngoscope Set, Miller, includes one each of Miller #0, #1, #2, #3, #4 blades; "AA" handle (#4066F) & "C" handle (#4065F) in hard case
#4079D	Fiberoptic Disposable Laryngoscope Set, Macintosh, includes two each of Macintosh #2, #3, #4 blades; "C" handle (#4065D) in soft case (#888)
#4089D	Fiberoptic Disposable Laryngoscope Set, Miller, includes two each of Miller #0, #1, #2, #3, blades; "C" handle (#4065D) in soft case (#888)

DISTRIBUTED BY:



AMERICAN DIAGNOSTIC CORPORATION

55 Commerce Drive
Hauppauge, New York 11788
tel: 1-800-232-2670 • fax: 631-273-9659

www.adctoday.com



Laryngcat118 rev 1
Printed in U.S.A.