


ADC® Digital Fingertip Pulse Oximeters

Diagnostix™ 2100 Digital Fingertip Pulse Oximeter

A portable, feature rich pulse oximeter that provides spot readings of SPO₂ and heart rate for the busy healthcare professional.

The 2100 features:

- Reads SPO₂ and Pulse Rate
- Accuracy: SPO₂ ±2% / HR ±5%BPM
- 6 changeable display modes
- 10 Display brightness settings
- Visual and audible alarms
- 25% smaller than competing models
- Bright 3 color OLED display
- Hinged battery cover is more secure
- Complete with Carry Lanyard and Black Rubberized Bumper (Optional orange bumper sold separately)
- 2 "AAA" Batteries
- Instruction Booklet
- 2 yr. Limited Warranty



Carry Lanyard



Rubberized Safety Bumper

Ordering Information

Item	Description
2100	Diagnostix™ Fingertip Pulse Oximeter
2100-1	Pulse Oximeter Carry Lanyard
2100-2BK	Pulse Oximeter Safety Bumper, Black
2100-2OR	Pulse Oximeter Safety Bumper, Orange




Advantage™ 2200 Digital Fingertip Pulse Oximeter

A portable finger tip pulse oximeter ideally suited for home, sport, or aviation use.

The 2200 features:

- Reads SPO₂ and Pulse Rate
- Accuracy: SPO₂ ±2% / HR ±5%BPM
- 2 changeable display modes
- Bright red LED display
- Carry Lanyard
- 2 "AAA" Batteries
- Instruction Booklet
- 1 yr. Limited Warranty



Carry Lanyard

Ordering Information

Item	Description
2200	Advantage™ Fingertip Pulse Oximeter
2100-1	Pulse Oximeter Carry Lanyard



1-800-ADC-2670
Visit us at www.adctoday.com

