


STETHOSCOPES

FEATURE CHART III



Adscope™ 600

The ultimate acoustic stethoscope

Features:

Benefits:

Stainless Steel Construction

Multi-frequency Design

Bi-lumen (two tube in one) design

ULTRA ADSOFT™ Threaded Eartips

Non-chill Diaphragm Rim

Binaural tubes fixed at 15°

Sculpted Chestpiece

Lifetime Warranty on Metal Parts

For the ultimate in acoustic performance and rugged durability

Allows detection of both Low and High Frequency sounds without the need to switch chestpiece sides

Greater stereo Separation without the risk of incidental rubbing

For the ultimate in acoustic performance and wearer comfort

For Patient comfort

For improved ergonomics

For a secure, stress free grip

Protects the owner's investment



ADSCOPE™ 600
Platinum Edition

Adscope™ 601

A convertible scope that combines the versatility of our Sprague
with the acoustic excellence of our 602

Features:

Benefits:

Stainless Steel Construction

Traditional Diaphragm/Bell Chestpiece

Convertible Chestpiece

Bi-lumen (two tube in one) design

ULTRA ADSOFT™ Threaded Eartips

Non-chill Diaphragm and Bell Rim

Binaural tubes fixed at 15°

Lifetime Warranty on Metal Parts

For the ultimate in acoustic performance and rugged durability

Time proven design for detecting the broadest possible range of high and low frequency sounds

Provides versatility for use with multiple patients

Greater stereo Separation without the risk of incidental rubbing

For the ultimate in acoustic performance and wearer comfort

For Patient comfort

For improved ergonomics

Protects the owner's investment



ADSCOPE™ 601
Convertible

Adscope™ 602

One of the finest acoustic scopes available, the 602 Cardiology model combines precision crafted performance with rugged durability

Features:

Benefits:

Stainless Steel Construction
Traditional Diaphragm/Bell Chestpiece

For the ultimate in acoustic performance and rugged durability

Time proven design for detecting the broadest possible range of high and low frequency sounds

Bi-lumen (two tube in one) design
ULTRA ADSOFT™ Threaded Eartips

Greater stereo Separation without the risk of incidental rubbing

For the ultimate in acoustic performance and wearer comfort

Non-chill Diaphragm Rim

For Patient comfort

Binaural tubes fixed at 15°

For improved ergonomics

Lifetime Warranty on Metal Parts

Protects the owner's investment



ADSCOPE™ 602
Cardiology

Adscope™ 603

One of our most popular acoustic scopes, the 603 combines outstanding acoustic performance with rugged durability

Features:

Benefits:

Stainless Steel Construction
Traditional Diaphragm/Bell Chestpiece

For the ultimate in acoustic performance and rugged durability

Time proven design for detecting the broadest possible range of high and low frequency sounds

Reinforcing Yoke in Binaural Assembly

Increases durability

ULTRA ADSOFT™ Threaded Eartips

For the ultimate in acoustic performance and wearer comfort

Non-chill Diaphragm and Bell Rim

For Patient comfort

Binaural tubes fixed at 15°

For improved ergonomics

Lifetime Warranty on Metal Parts

Protects the owner's investment

Wide selection of colors

For personalization, or departmental coding



ADSCOPE™ 603
Acoustic

Adscope™ 604 / 605

Legendary acoustics and performance in our finest Adscope™ in scaled down versions for exceptional functionality for pediatric and infant care

Features:

Benefits:

Stainless Steel Construction

Traditional Diaphragm/Bell Chestpiece

Reinforcing Yoke in Binaural Assembly

ULTRA ADSOFT™ Threaded Eartips

Non-chill Diaphragm and Bell Rim

Binaural tubes fixed at 15°

Lifetime Warranty on Metal Parts

Wide selection of colors

For the ultimate in acoustic performance and rugged durability

Time proven design for detecting the broadest possible range of high and low frequency sounds

Increases durability

For the ultimate in acoustic performance and wearer comfort

For Patient comfort

For improved ergonomics

Protects the owner's investment

For personalization, or departmental coding



ADSCOPE™ 604 / 605
Acoustic Pediatric / Infant

Adscope-Lite™ 609

The Adscope-Lite™ is the only affordably priced stethoscope to combine superior acoustic response with an ultra lightweight design

Features:

Benefits:

Aluminum Alloy Construction

Traditional Diaphragm/Bell Chestpiece

Reinforcing Yoke in Binaural Assembly

ULTRA ADSOFT™ Threaded Eartips

Non-chill Diaphragm and Bell Rim

Binaural tubes fixed at 15°

Wide selection of colors

Balances acoustic performance with lightweight comfort

Time proven design for detecting the broadest possible range of high and low frequency sounds

Increases durability

For the ultimate in acoustic performance and wearer comfort

For Patient comfort

For improved ergonomics

For personalization, or departmental coding



ADSCOPE-LITE™ 609

Adscope™ 612 / 614

The lightweight version of our world famous multi-frequency scopes
now in both adult and pediatric sizes

Features:

Benefits:

Zinc Construction

Provides the acoustic performance of Stainless steel without the cost

Multi-frequency Design

Allows detection of both Low and High Frequency sounds without the need to switch chestpiece sides

Reinforcing Yoke in Binaural Assembly

Increases durability

ULTRA ADSOFT™ Threaded Eartips

For the ultimate in acoustic performance and wearer comfort

Non-chill Diaphragm Rim

For Patient comfort

Binaural tubes fixed at 15°

For improved ergonomics

Sculpted Chestpiece

For a secure, stress free grip

Lifetime Warranty on Metal Parts

Protects the owner's investment



ADSCOPE™ 612 / 614
Platinum Edition Adscope-Lite™

Adscope™ 615

A convertible scope that combines the versatility of our Sprague
with the acoustic excellence of our 602

Features:

Benefits:

Stainless Steel Construction

For the ultimate in acoustic performance and rugged durability

Traditional Diaphragm/Bell Chestpiece

Time proven design for detecting the broadest possible range of high and low frequency sounds

Convertible Chestpiece

Provides versatility for use with multiple patients

Bi-lumen (two tube in one) design

Greater stereo Separation without the risk of incidental rubbing

ULTRA ADSOFT™ Threaded Eartips

For the ultimate in acoustic performance and wearer comfort

Non-chill Diaphragm Rim

For Patient comfort

Binaural tubes fixed at 15°

For improved ergonomics

Lifetime Warranty on Metal Parts

Protects the owner's investment



ADSCOPE™ 615
Platinum Professional Edition

Adscope™ 617 / 618

Our newest multi-frequency scope with a fresh face!

Features:

Benefits:

Construction

Distinctive molded resin animal face chestpiece in choice of adult or pediatric sizes

Multi-frequency Design

Multi-frequency design combines the acoustic response of a traditional bell and diaphragm into a single sided chestpiece

Reinforcing Yoke in Binaural Assembly

Increases durability

ULTRA ADSOFT™ Threaded Eartips

For the ultimate in acoustic performance and wearer comfort

Gasketed Diaphragm

Switches between bell-like low frequency and diaphragm-like high frequency performance simply by altering pressure

Aluminum Binaural tubes fixed at 15°

For improved ergonomics

Snap-on Chestpiece

7 interchangeable Animal face snap-ons include: Panda, Frog, Koala, Deer, Monkey, Tiger, and Bear

Five Year Warranty

Protects the owner's investment



ADSCOPE™ 617 / 618
Stethoscopes

Adscope™ 640BK / 641

Our exclusive design features an improved valve mechanism that virtually eliminates acoustic leakage and a chestpiece key to facilitate after market service

Features:	Benefits:
Zinc chestpiece Construction	Provides the acoustic performance of Stainless steel without the cost
Convertible Chestpiece	Provides versatility for use with multiple patients
Traditional Diaphragm/Bell Chestpiece	Time proven design for detecting the broadest possible range of high and low frequency sounds
Double Tube configuration	Provides enhanced acoustic performance and stereo separation
Patented Valve mechanism	Virtually eliminates acoustic leakage found on competing sprague designs
ULTRA ADSOFT™ Threaded Eartips	For the ultimate in acoustic performance and wearer comfort (641 only)
Gasketed drum	Improves acoustic performance and security of chestpiece attachments
Chestpiece key	Facilitates removal of worn or damaged diaphragms.
Wide selection of colors	For personalization, or departmental coding (641 only)



ADSCOPE™ 640BK / 641
Sprague Rappaport

Adscope™ 647 Sprague-1

ADC's legendary Sprague with a lightweight single tube design for greater comfort and convenience.

Features:	Benefits:
Zinc chestpiece Construction	Provides the acoustic performance of Stainless steel without the cost
Convertible Chestpiece	Provides versatility for use with multiple patients
Traditional Diaphragm/Bell Chestpiece	Time proven design for detecting the broadest possible range of high and low frequency sounds
40% thicker Y Tubing	Reduces the weight and bulk of traditional spragues
ULTRA ADSOFT™ Threaded Eartips	For the ultimate in acoustic performance and wearer comfort
Patented Valve mechanism	Virtually eliminates acoustic leakage found on competing sprague designs
Gasketed drum	Improves acoustic performance and security of chestpiece attachments
Chestpiece key	Facilitates removal of worn or damaged diaphragms.



ADSCOPE™ 647
Sprague-1™

Adscope™ 657

The next generation in electronic stethoscopy, combining impressive acoustic amplification, electronic noise filtering and comfortable ergonomic design

Features:	Benefits:
3 listening modes	For the broadest possible frequency response
Contact microphone	ONLY picks up sounds detected from what the chestpiece touches eliminating a lot of ambient noise
8 step volume control	Allows very precise control of desired volume
16x amplification	Provides 24db boost over traditional acoustic scopes
Remembers last mode and volume	Practitioner need not reset the volume control or mode each time the instrument is used
"AAA" Alkaline Battery	Provides up to 180 hours continuous use
3 minute Auto Off	Conserves battery life
Auto Off Cancel	Touching ANY button resets the auto off counter to 0
Low battery warning LED	Light changes from amber to red to indicate weakening battery
Binaural tubes fixed at 15°	For improved ergonomics
ULTRA ADSOFT™ Threaded Eartips	For the ultimate in acoustic performance and wearer comfort



ADSCOPE™ 657
Electronic Scope

Adscope™ 655 VistaScope

A unique chestpiece design with 8 interchangeable designer diaphragms which provides up to twice the amplification of ordinary acoustic scopes

Features:	Benefits:
Oversized Acrylic Chestpiece	Improves acoustic performance by creating a resonant chamber
Interchangeable diaphragms which show through chestpiece	Allow the practitioner to personalize their scope
Reinforcing Yoke in Binaural Assembly	Increases durability
German Stainless Steel binaural	Rugged durability without the extra weight
ULTRA ADSOFT™ Threaded Eartips	For the ultimate in acoustic performance and wearer comfort
Binaural tubes fixed at 15°	For improved ergonomics
Lifetime Warranty on Metal Parts	Protects the owner's investment



ADSCOPE™ 655
VistaScope

Proscope™ 660

ADC Nurse style Proscope™ combines lightweight comfort, sensitive acoustics, and handsome styling, with truly affordable pricing

Features:

Benefits:

Lightweight aluminum construction

Durable, yet extraordinarily comfortable, weighs just 3 oz.

**Extra long binaural
(up to 1/2" longer than competing brands)**

Adjustable providing greater comfort for the practitioner.
Overall length 32"

Proprietary diaphragm/rim design

Enhances acoustic performance by up to 50%.
Ideally suited for blood pressure measurement

Non-chill diaphragm Rim


For Patient comfort

Spare Mushroom Eartips

For practitioner wearing comfort

Wide selection of colors

Choice of 21 colors for personalization, or departmental coding



PROSCOPE™ 660
Nurse Style

Proscope™ 670

ADC Proscope™ combines lightweight comfort, sensitive acoustics, in a versatile Dual-Head scope

Features:

Benefits:

Combination Diaphragm/Bell Chestpiece

Time proven design for detecting the broadest possible range of high and low frequency sounds

Lightweight aluminum construction

Durable, yet extraordinarily comfortable

**Extra long binaural
(up to 1/2" longer than competing brands)**

Provides greater comfort for the practitioner
Overall length 32"

Proprietary diaphragm/rim design

Enhances acoustic performance by up to 50%

Non-chill diaphragm Rim


For Patient comfort

Spare Mushroom Eartips

For practitioner wearing comfort

Wide selection of colors

For personalization, or departmental coding



PROSCOPE™ 670
Dual-Head

Proscope™ 662

Specialty institutional traditional Bowles Scope

Features:

Benefits:

Lightweight aluminum binaural

Durable, yet lighter weight than brass

**Extra long binaural
(up to 1/2" longer than competing brands)**

Provides greater comfort for the practitioner

Non-chill diaphragm Rim

For Patient comfort

Spare Mushroom Eartips

For practitioner wearing comfort



PROSCOPE™ 662
Bowles Scope

Proscopes™ 664 / 665

The first disposable scopes with performance so good you might want to keep it

PROSCOPE™ 664
Disposable



Features:

Benefits:

Lightweight aluminum binaural (664 Only)

Durable, yet lighter weight than brass

Completely disposable stethoscope with no metal parts (665 Only)

Easily incinerated

Ultra lightweight

For comfort and cost effectiveness

Multi-level packaging (ea in poly bag, 10/ctn, 10 ctn/case)

Improved stock control

**664 available in Yellow
665 available in Yellow and Royal Blue**

For infection control areas



PROSCOPE™ 665
Disposable

Proscope™ 675

Institutional Pediatric Scope

Features:

Benefits:

Traditional Diaphragm/Bell Chestpiece

Time proven design for detecting the broadest possible range of high and low frequency sounds

Specially sized chestpiece

Ideally suited for pediatric applications

Lightweight aluminum construction THROUGHOUT

Durable, yet extraordinarily comfortable

Extra long binaural (up to 1/2" longer than competing brands)

Provides greater comfort for the practitioner

Non-chill diaphragm Rim

For Patient comfort

Spare Mushroom Eartips

For practitioner wearing comfort

Wide selection of colors

For personalization, or departmental coding



PROSCOPE™ 675
Institutional Pediatric Scope

Proscope™ 676

Institutional Infant Scope

Features:

Benefits:

Traditional Diaphragm/Bell Chestpiece

Time proven design for detecting the broadest possible range of high and low frequency sounds

Specially sized chestpiece

Ideally suited for infant applications

Lightweight aluminum binaural

Durable, yet lighter weight than brass

Extra long binaural (up to 1/2" longer than competing brands)

Provides greater comfort for the practitioner

Non-chill Diaphragm Rim

For Patient comfort

Spare Mushroom Eartips

For practitioner wearing comfort

Wide selection of colors

For personalization, or departmental coding



PROSCOPE™ 676
Institutional Infant Scope

ADC® Model / Brand		600 Adscope™	601 Adscope™	602 Adscope™	603 Adscope™	604 Adscope™	605 Adscope™	609 Adscope™	612 Adscope™	614 Adscope™	615 Adscope™	617 Adscope™ • 618 Adscope™
FEATURES	Product Name	Platinum Multi-Frequency	Convertible Cardiology	Cardiology	Adult Stainless	Pediatric Stainless	Infant Stainless	Adult Lightweight	Platinum-Lite Adult	Platinum-lite Pediatric	Platinum Professional	Multi-Frequency
	Type	Multit-frequency	Cardiology	Cardiology	Traditional	Traditional	Traditional	Traditional	Multit-frequency	Multit-frequency	Multit-frequency	Multit-frequency
	Chestpiece Size	Adult	Adult/Ped	Adult	Adult	Pediatric	Infant	Adult	Adult	Pediatric	Adult	Adult • Pediatric
	Diaphragm Dim	1 3/4"	1 3/4"	1 3/4"	1 3/4"	1 3/8"	1 1/8"	1 3/4"	1 1/2"	1"	1 7/8" x 2 3/16"	1 7/8" • 1 1/2"
	Bell Dim	NA	1 1/2"	1 3/8"	1 3/8"	1 1/8"	7/8"	1 3/8"	NA	NA	NA	NA
	Chestpiece type	Multit-frequency	Diaphragm/Bell	Diaphragm/Bell	Diaphragm/Bell	Diaphragm/Bell	Diaphragm/Bell	Diaphragm/Bell	Multi-frequency	Multi-frequency	Multi-frequency	Multi-frequency
	Binaural Spring type	Internal	Internal	Internal	Internal	Internal	Internal	Internal	Internal	Internal	Internal	Internal
	Tubing Type	Bi-lumen	Bi-lumen	Bi-lumen	Single lumen	Single lumen	Single lumen	Single lumen	Single lumen	Single lumen	Single lumen	Single lumen
	Chestpiece Material	Stainless	Stainless	Stainless	Stainless	Stainless	Stainless	Aluminum	Zinc	Zinc	Zinc	Molded Resin
	Binaural Material	Stainless	Stainless	Stainless	Stainless	Stainless	Stainless	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum
	Tubing Length	19"	19"	16"	22"	22"	22"	22"	22"	22"	22"	22"
	Overall Length	27"	28"	25.5"	31"	30.5"	30.5"	31"	30.5"	30.5"	30.5"	30"
	Weight	7.5 oz	9oz (A), 10oz (P)	7.5 oz	6.5 oz	5.25 oz	4.5 oz	4.5 oz	5 oz	3.9 oz	5.4 oz	4 oz • 3 oz
	Spare Eartips	602-07,602-09	602-07,602-09	602-07,602-09	602-07,602-09	602-07,602-09	602-07,602-09	602-09	602-07, 602-09	602-07, 602-08	602-07, 602-08	NA
	Spare Diaphragm	601-03	601-03	602-03	602-03	604-03	605-03	609-03	601-03	614-03	615-03	600-03NG, 614-03
INTENDED USES	Scope ID tag	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
	Tubing Colors	4	4	3	13	4	4	18	4	4	7	3
	Warranty	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	5 Year	Lifetime	Lifetime	Lifetime	5 Year
	Cardiology	•	•	•								
	Respiratory Specialist	•	•	•								
	General Practitioner	•	•	•	•						•	•
	ER	•	•	•	•						•	
	Medical Student	•	•	•	•						•	
	Pediatric Specialitsts					•	•			•		•
	Nurse				•			•	•		•	•
	Critical Care Nurses	•	•	•	•						•	
	Entry Level Nurses							•	•			•
	EMT	•	•	•	•			•	•		•	•
	Institutional	•	•									

ADC® Model / Brand		646 Adscope™	647 Adscope™	655 Adscope™	657 Adscope™	660 Proscope™	662 Proscope™	664 Proscope™	665 Proscope™	670 Proscope™	675 Proscope™	676 Proscope™
FEATURES	Product Name	Sprague	Sprague	Vista Scope	Adscope 657	Nurse	Bowles	SPU	SPU	Dual Head	Dual Head Pediatric	Dual Head Infant
	Type	Sprague	Sprague	Specialty	Electronic	Institutional	Institutional	Institutional	Institutional	Institutional	Institutional	Institutional
	Chestpiece Size	Multi	Multi	Adult	Multi	Adult	Adult	Adult	Adult	Adult	Pediatric	Infant
	Diaphragm Dim	1 1/4", 1	1 1/4", 1	1 7/8"	1 3/4"	1 3/4"	1 5/8"	1 3/4"	1 3/4"	1 3/4"	1 1/4"	.98" (25mm)
	Bell Dim	1 1/4", 1", 3/4"	1 1/4", 1", 3/4"	NA	NA	NA	NA	NA	NA	1 1/8"	3/4"	.98" (25mm)
	Chestpiece type	Diaphragm/Bell	Diaphragm/Bell	Diaphragm	Electronic	Diaphragm	Diaphragm	Diaphragm	Diaphragm	Diaphragm/Bell	Diaphragm/Bell	Diaphragm/Bell
	Binaural Spring type	External	External	Internal	Internal	External	External	External	External	External	External	External
	Tubing Type	Two Tube	Two Tube	Single Lumen	Single Lumen	Single Lumen Y	Single Lumen Y	Single Lumen Y	Single Lumen Y	Single Lumen Y	Single Lumen Y	Single Lumen Y
	Chestpiece Material	Zinc	Zinc	Acrylic	Various Materials	Aluminum	Zinc	PVC	PVC	Aluminum	Aluminum	Aluminum
	Binaural Material	Brass	Brass	Stainless	Aluminum	Aluminum	Aluminum	Aluminum	ABS	Aluminum	Aluminum	Aluminum
	Tubing Length	22"	22"	22"	22"	22"	22"	22"	22"	22"	22"	22"
	Overall Length	30"	30"	31.5"	28"	32"	32.5"	32"	32"	32.5"	31.5"	32"
	Weight	8 oz	6.2 oz	5 oz	6 oz	3 oz	3.25 oz	2.5 oz	1.5 oz	3.5 oz	3.5 oz	4 oz
	Spare Eartips	640-10, 600-06BK	640-10, 600-06BK	602-07, 602-09	602-07	600-06BK	600-06BK	600-06BK	600-06BK	600-06BK	600-06BK	600-06BK
	Spare Diaphragm	640-03 640-03S	640-03 640-03S	655-03K	NA	NA	NA	NA	NA	NA	NA	NA
INTENDED USES	Scope ID tag	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO
	Tubing Colors	1	5	6	1	21	1	1	2	21	21	21
	Warranty	5 Year	5 Year	Lifetime	2 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year
	Cardiology				•							
	Respiratory Specialist				•							
	General Practitioner				•							
	ER				•							
	Medical Student				•							
	Pediatric Specialitsts											
	Nurse	•	•	•								
	Critical Care Nurses				•							
	Entry Level Nurses	•	•	•						•		
	EMT	•	•	•	•					•		
	Institutional	•	•	•		•	•	•	•	•	•	•

STETHOSCOPE COLOR CHART

	BK	BB	BD	CA	DG	FG	FI	FKW	FL	FM	FPB	FPL	FRB	FS	FV	G	GD	GR	HP	LB	LV	M	N	NB	NGR	NO	NP	NY	O	OD	P	R	RB	SF	ST	TL	TQ	V	W	Y		
Adscope™ 600																																										
Adscope™ 600GP																																										
Adscope™ 601																																										
Adscope™ 602																																										
Adscope™ 603																																										
Adscope™ 604																																										
Adscope™ 605																																										
Adscope™ 609																																										
Adscope™ 612																																										
Adscope™ 614																																										
Adscope™ 615																																										
Adimals™ 617/618																																										
Proscope™ 640																																										
Adscope™ 641																																										
Adscope™ 645																																										
Adscope™ 646ST																																										
Adscope™ 647																																										
Adscope™ 655																																										
Adscope™ 657BK																																										
Proscope™ 660																																										
Proscope™ 661																																										
Proscope™ 662																																										
Proscope™ 664																																										
Proscope™ 665																																										
Proscope™ 670																																										
Proscope™ 671																																										
Proscope™ 675																																										
Proscope™ 676																																										
Scope Tag 697BK																																										
Scope Tag 698BK																																										

*Colors shown may vary from actual color.

STETHOSCOPE CROSS REFERENCE CHART

ADC	Description	Baum	BV	Grafco	Littmann	Mabis	Medline	Prestige	Trimline	Welch Allyn
600	Platinum Multifrequency		13-450	04-409	2159, 216X	10-410-020				
600GP	Platinum Multifrequency GP							126-G		
601	Convertible Cardiology				312X			112		5079-32X
602	Cardiology	273	13-74X			10-102-0X0	MDS9250	128, 129	59101	5079-1XX
603	Adult Stainless	272	13-70X, 401		220X	10-404-0X0	MDS922X	126, 127	59105	5079-13X
604	Pediatric Stainless		13-71X		2113	10-406-0X5	MDS9220		59111	5079-14X
605	Infant Stainless		13-72X		2114	10-408-0X3	MDS921X		59115	5079-22X
609	Adult Lightweight				2450			121	59103	
612	Platinum-Lite, Adult		13-250		2290	10-412-020				
614	Platinum-Lite, Pediatric									
615	Platinum Pro		13-410		4470					
640	Sprague			04-600						
641	Sprague	270	416-22			10-4140X0	MDS92630X	122, S132		5079-28
645	Sprague, GP							122-G		
646	Sprague, Stealth									
647	Sprague, One Tube							S124, S134		
655	VistaScope		13-320, 510			10-460-3X0		S107		
657	Adscope™ 657				3000, 4100WS	10-400-020		130	59001	
660	Proscope™, Nurse		411	04-300		10-424-0X0, 10-428-0X0	MDS92610X	106	5925X	
662	Proscope™, Bowles		407			10-440-020	MDS9515			
664	Proscope™, SPU		13-176H							
665	Proscope™, SPU		300	04-722		10-448-130	MDSG543	100	59801	1746X
670	Proscope™, Dual Head		412	04-400		10-422-0X0, 10-426-0X0	MDS92620X	108, S125	5920X	
671	Proscope™, Teaching		405	04-540		10-446-020	MDS9750	S122-T		
675	Proscope™, Dual Head, Ped		13-300	04-507		10-432-0X5	MDS9557	S108-P		
676	Proscope™, Dual Head, Infant		13-310	04-512		10-434-0X2	MDS9565	S108-I		



ADC
55 Commerce Drive
Hauppauge, NY 11788

tel: 631-273-9600
toll free: 1-800-232-2670
fax: 631-273-9659

www.adctoday.com
info@adctoday.com



ADC (UK) Ltd.
Unit 6, PO14 1TH
United Kingdom



6 34782 05979 8