


DIGITAL THERMOMETERS

FEATURE CHART III



Thermometer Products Feature Chart



ADC Item #	412 Adtemp™ I	413B/R Adtemp™ II	415/415FL Adtemp™ IV	418 Adtemp™ V	419 Adtemp™ Hypothermia	421 Adtemp™ Tympanic	422 Adtemp™ Veterinary	423 Adtemp™ VII	424 Adtemp™ Tympanic	427 Adtemp™ TempleTouch
Type	Oral • Rectal Axillary	Oral • Rectal Axillary	Oral • Rectal Axillary	Oral • Rectal Axillary	Oral • Rectal Axillary	Ear	Rectal	Oral • Rectal Axillary	Ear	Temple
Scale	Dual (°F/°C)	Dual (°F/°C)	Dual (°F/°C)	Dual (°F/°C)	Dual (°F/°C)	Dual (°F/°C)	Dual (°F/°C)	Dual (°F/°C)	Dual (°F/°C)	Dual (°F/°C)
Range °F	90.0-109.9	90.0-109.9	90.0-109.9	90.0-109.9	78.0-111.9	89.6-108.0	90.0-111.9	90.0-109.9	86.6-108.0	95.0-107.6
Range °C	32.2-43.9	32.2-43.9	32.2-43.9	32.2-43.9	25.5-44.3	32-42.2	32.0-43.9	32.0-43.9	32.0-42.2	35.0-42.0
Accuracy °F	0.2	0.2	0.2	0.2	0.2	0.4	0.2	0.2	0.2	0.2
Accuracy °C	0.1	0.1	0.1	0.1	0.1	0.2	0.1	0.1	0.1	0.1
Response Time	60 sec	60 sec	10/18 sec	8 sec	10/18 sec	1 sec	10 sec	60 sec	1 sec	6 sec
Beeper	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Flex Probe	No	No	No / Yes	Yes	No	No	No	No	No	No
Need Sheaths	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes	No
Warranty	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year
Battery Type	1.55v (LR41)	1.55v (LR41)	1.55v (LR41)	1.55v (LR41)	1.55v (LR41)	3 v	1.55v (LR41)	1.55v (LR41)	1.55v (CR2032)	2 x "AAA"
Latex Content	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free

Thermometer Products Cross Reference Chart

Omron	MC-110	MC-130MO ³	MC-106	MC-206	N/A	N/A	N/A	N/A	MC-514	N/A
	N/A	MC-130MR ⁴	N/A	N/A	N/A	N/A	N/A	N/A	MC-505	N/A
Prestige	N/A	DT	N/A	DT-3	N/A	N/A	N/A	N/A	DT-23	N/A
	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	N/A	DT-6 ⁵	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Mabis	15-689-000	15-719-000 ³	15-720-000	15-736-000	N/A	18-207-000	N/A	15-600-000 ¹	18-200-000	N/A
	N/A	15-749-000 ⁴	15-732-000	N/A	N/A	N/A	N/A	15-604-000 ²	N/A	N/A
	N/A	15-713-000	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	N/A	15-691-000 ⁹	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Lumiscope	N/A	2013	2210	N/A	N/A	2215	N/A	N/A	N/A	N/A
BD	N/A	524100	N/A	524101	N/A	N/A	N/A	N/A	524-972	N/A
Vicks (Kaz)	820	N/A	V911	V965	N/A	N/A	N/A	V900	V971N	V977
	821-24	N/A	V912	V966	N/A	N/A	N/A	V901	N/A	N/A
	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Timex	N/A	N/A	N/A	80002 ¹⁰	N/A	N/A	N/A	N/A	80010	N/A
	N/A	N/A	N/A	80006 ¹⁰	N/A	N/A	N/A	N/A	N/A	N/A
LifeSource	N/A	DT-703	DT-704	DT-705	N/A	UT-302	N/A	N/A	N/A	N/A
Braun	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	IRT-3020CO	N/A

Notes:

- 1) Fahrenheit only
2) Celsius only

- 5) Available in 6 colors
9) Not waterproof

- 10) Response times differ from ADC,
also has back light and curved probe tip

Adtemp™ I Digital Thermometer #412



Features:	Benefits:
Dual Scale	For readability in Fahrenheit or Celsius
Water Proof Probe	Easy cleaning
Last Reading Memory	Aids in comparisons amongst readings
2 Year Warranty	Consumer protection
Latex-free	No risk of allergic reactions
Bulk packaging	Creates low cost
Kitted packaging	Complete with 20 sheaths
Compact Size	Comfortable use
Auto-off	Conserves battery life
1.55v Battery	Provides up to 1,500 measurements
Integral carry case	Keeps instrument sanitary

Sold as:

#412-00	Bulk Packaging No sheaths, polybag, no label on units 20/Cartron, 240/Case
#412K	Kitted Packaging 20 sheaths, 2 pocket vinyl wallet 1/Box, 12/Cartron, 144 Case

Adtemp™ II Digital Thermometer #413B or #413R



Features:	Benefits:
Dual Scale	For readability in Fahrenheit or Celsius
Water Proof Probe	Easy cleaning
Last Reading Memory	Aids in comparisons amongst readings
2 Year Warranty	Consumer protection
Latex-free	No risk of allergic reactions
Gift box packaging	Complete with 5 sheaths, displays product
Bulk packaging	Creates low cost
Kitted packaging	Complete with 20 sheaths
Compact Size	Comfortable use
Auto-off	Conserves battery life
1.55v Battery	Provides up to 1,500 measurements
Integral carry case	Keeps instrument sanitary
Color coded blue or red	Referencing oral or rectal use, same functionality
Beeps when complete	Lets user know when finished reading temperature
Replaceable battery	Encourages longer use of instrument

Sold as:

#413B	Oral, Blue Cap, Gift Box Packaging 5 sheaths 1/Box, 12/Cartron, 240/Case
#413B-00	Oral, Blue Cap, Bulk Packaging 5 sheaths, polybag, no label on units 25/Cartron, 500/Case
#413BK	Oral, Blue Cap, Kitted Packaging 20 sheaths, 2 pocket vinyl wallet 1/Box, 12/Cartron, 144/Case
#413R	Rectal, Red Cap, Gift Box Packaging 5 sheaths 1/Box, 12/Cartron, 240/Case
#413R-00	Rectal, Red Cap, Bulk Packaging 5 sheaths, polybag, no label on units 25/Cartron, 500/Case
#413RK	Rectal, Red Cap, Kitted Packaging 20 sheaths, 2 pocket vinyl wallet 1/Box, 12/Cartron, 144/Case

Adtemp™ IV Digital Thermometer #415 / #415FL

Dual Scale	For readability in Fahrenheit or Celsius
Water Proof Probe	Easy cleaning
Last Reading Memory	Aids in comparisons amongst readings
Jumbo display	Easy reading
2 Year Warranty	Consumer protection
Latex-free	No risk of allergic reactions
Flexible tip	Patient safety (415FL model only)
Gift box packaging	Complete with 5 sheaths, displays product
Kitted packaging	Complete with 20 sheaths
Auto-off	Conserves battery life
1.55v Battery	Provides up to 1,500 measurements
Integral carry case	Keeps instrument sanitary
Beeps when complete	Lets user know when finished reading temperature
Replaceable battery	Encourages longer use of instrument
10 second fast read	Earlier response time than standard model



Adtemp™ IV 415



Adtemp™ IV 415FL

Sold as:

#415	Gift Box Packaging 5 sheaths 1/Box, 6/Case, 72/Case
#415FL	Gift Box Packaging 5 sheaths 1/Box, 6/Case, 72/Case
#415K	Kitted Packaging 20 sheaths, 2 pocket vinyl wallet 1/Box, 12/Case, 144/Case
#415FLK	Kitted Packaging 20 sheaths, 2 pocket vinyl wallet 1/Box, 12/Case, 144/Case

Adtemp™ V Digital Thermometer #418

Dual Scale	For readability in Fahrenheit or Celsius
Water Proof Probe	Easy cleaning
Last Reading Memory	Aids in comparisons amongst readings
Jumbo display	Easy reading
2 Year Warranty	Consumer protection
Latex-free	No risk of allergic reactions
Gift box packaging	Complete with 5 sheaths, displays product
Auto-off	Conserves battery life
1.55v Battery	Provides up to 1,500 measurements
Integral carry case	Keeps instrument sanitary
Beeps when complete	Lets user know when finished reading temperature
Replaceable battery	Encourages longer use of instrument
Flexible tip	Patient safety
8 second fast read	Earlier response time than standard model



Sold as:

#418	Gift Box Packaging 5 sheaths 1/Box, 12/Case, 240/Case
------	---

Adtemp™ Hypothermia Digital Thermometer #419

Features:	Benefits:
Extended range	Detects severe low temperatures
Auto-off	Conserves battery life
First digital	Mercury free
Replaceable battery	Encourages longer use of instrument
Integral carry case	Keeps instrument sanitary
Bulk packaging	Creates low cost
Dual Scale	For readability in Fahrenheit or Celsius
Water Proof Probe	Easy cleaning
Last Reading Memory	Aids in comparisons amongst readings
Jumbo display	Easy reading
2 Year Warranty	Consumer protection
Latex-free	No risk of allergic reactions
1.55v Battery	Provides up to 1,500 measurements
Beeps when complete	Lets user know when finished reading temperature




Sold as:

#419 Bulk Packaging
1/Polybag, 25/Cartron, 500/Case

Adtemp™ Tympanic Digital Thermometer #421

Features:	Benefits:
Measured in ear	Patient comfort
No probe covers needed	Saves on additional costs and more convenient
Dual Scale	For readability in Fahrenheit or Celsius
Water Proof Probe	Easy cleaning
12X Reading Memory	Aids in comparisons amongst readings
Jumbo display	Easy reading
2 Year Warranty	Consumer protection
Latex-free	No risk of allergic reactions
Gift Box packaging	Displays product
Auto-off	Conserves battery life
Beeps when complete	Lets user know when finished reading temperature
Replaceable battery	Encourages longer use of instrument
1 second reading	Fastest response time
Integral carry case	Easy for storage and travel
Ultra Compact Size	Comfortable use



Sold as:

#421 Gift Box Packaging
1/Box, 6/Cartron, 48/Case

Adtemp™ Veterinary Digital Thermometer #422

Features:	Benefits:
Veterinary applications	Patient comfort
Specially designed housing with lanyard attachment	Enhanced removal
Dual Scale	For readability in Fahrenheit or Celsius
Water Proof Probe	Easy cleaning
Last Reading Memory	Aids in comparisons amongst readings
Jumbo display	Easy reading
2 Year Warranty	Consumer protection
Latex-free	No risk of allergic reactions
Bulk packaging	Creates low cost
Auto-off	Conserves battery life
Beeps when complete	Lets user know when finished reading temperature
1.55v battery	Provides up to 1,500 measurements
Integral carry case	Keeps instrument sanitary
Replaceable battery	Encourages longer use of instrument
10 second fast read	Earlier response time than standard model



Sold as:

#422 Bulk Packaging
1/Polybag, 25/Case, 500/Case

Adtemp™ VII Digital Thermometer #423

Features:	Benefits:
Dual Scale	For readability in Fahrenheit or Celsius
Water Proof Probe	Easy cleaning
Last Reading Memory	Aids in comparisons amongst readings
Jumbo display	Easy reading
2 Year Warranty	Consumer protection
Latex-free	No risk of allergic reactions
Gift Box Packaging	Complete with 5 sheaths, displays product
Kitted packaging	Complete with 20 sheaths
Auto-off	Conserves battery life
1.55v Battery	Provides up to 1,500 measurements
Integral carry case	Keeps instrument sanitary
Beeps when complete	Lets user know when finished reading temperature
Replaceable battery	Encourages longer use of instrument



Sold as:

#423 Gift Box Packaging
5 Sheaths
1/Box, 24/Case, 288/Case

#423K Kitted Packaging
20 sheaths, 2 pocket vinyl wallet
1/Box, 12/Case, 144 Case

Adtemp™ Tympanic Digital Thermometer #424

Features:	Benefits:
Measured in ear	Patient comfort
Includes probe covers	Reduces the risk of cross contamination
Dual Scale	For readability in Fahrenheit or Celsius
Water Proof Probe	Easy cleaning
Last Reading Memory	Aids in comparisons amongst readings
Jumbo display	Easy reading
2 Year Warranty	Consumer protection
Latex-free	No risk of allergic reactions
Gift Box packaging	Protects product
Auto-off	Conserves battery life
Beeps when complete	Lets user know when finished reading temperature
Replaceable battery	Encourages longer use of instrument
1 second reading	fastest response time
Desktop storage caddy	Easy for storage and travel
Compact ergonomic size	Comfortable use



Sold as:

#424	Gift Box Packaging 1/Box, 6/Cartron, 48/Case
#425-40	Additional Probe Cover Dispenser (40 Covers per Dispenser) 1/Box, 24/Cartron, 96/Case

Adtemp™ Temple Touch Thermometer #427

Features:	Benefits:
Measured at the temple	For patient comfort
No probe covers needed	Saves on additional costs and more convenient
Dual Scale	For readability in Fahrenheit or Celsius
Water Proof Probe	Easy cleaning
Last Reading Memory	Aids in comparisons amongst readings
Jumbo display	Easy reading
2 Year Warranty	Consumer protection
Latex-free	No risk of allergic reactions
Gift Box Packaging	Protects product
Clamshell Packaging	Displays product
Auto-off	Conserves battery life
Beeps when complete	Lets user know when finished reading temperature
Replaceable Batteries	Encourages longer use of instrument
6 second reading	Fastest response time
Backlit display	Easy reading
Ultra compact size	Comfortable use



Sold as:

#427	Gift Box Packaging 1/Box, 12/Cartron, 144/Case
#427Q	Display Packaging 1/Box, 12/Cartron, 36/Case



ADC
55 Commerce Drive
Hauppauge, NY 11788



ADC (UK) Ltd.
Unit 6, PO14 1TH
United Kingdom

tel: 631-273-9600
toll free: 1-800-232-2670
fax: 631-273-9659

www.adctoday.com
info@adctoday.com

