

# Adtemp™ 424 Tympanic Thermometer





Large Backlit LCD Display



Audible "fever" alarm for temperatures above 99.5°F/37.5°C.



Switchable Dual scale Measures in either °F or °C modes.



Compact Ergonomic Design & Storage Caddy



Convenient Probe Cover Dispenser Contains 20 Covers

## ADTEMP™ 424 INFRARED EAR THERMOMETER

Easy to use tympanic thermometer measures core temperature instantly. Convenient, economical probe covers make it the ideal instrument for use on all types of patients.

### ADTEMP™ 424 Tympanic thermometer features:

- 1 second measurement using scan/peak technology
- Large 6mm illuminated display
- Temp range from 89.6°F - 108°F +/- .2°F (32°C - 42.2°C +/- .1°C)
- F/C switchable display
- Fever alarm at 99.5°F (37.5°C)
- Last reading memory
- Audible beep when measurement complete
- Auto off conserves battery strength
- Compact, ergonomic design - measures 5½" (13.97cm) long, weighs 1.6oz. (45.3g) including battery
- Powered by 1.55 (CR20320) Lithium battery (included)
- Includes 20 probe covers in convenient dispenser
- Desktop storage caddy



ORDERING INFORMATION	
#424	Adtemp™ Infrared Ear Thermometer
#425-40	Adtemp™ Tympanic Sheaths, 40/pkg



1-800-ADC-2670  
Visit us at [www.adctoday.com](http://www.adctoday.com)



# Adtemp™ 427 Temple Touch Thermometer

TM




## ADTEMP™ 427 TEMPLE TOUCH THERMOMETER

Provides an accurate 6 second non-invasive temperature reading at the forehead making it ideal for pediatric patients. Cost effective design makes it ideal for institutional use.

### ADTEMP™ 427 Temple Touch thermometer features:

- Non-invasive forehead temperature measurement using patented conductive technology delivers clinically accurate readings in about 6 seconds
- Reduces the risk of cross contamination and avoids the need for costly probe covers
- Ideal for the hospital, clinic, or physician office
- Auto off function conserves battery life
- Audible tone when measurement is complete
- Last reading memory
- Dual Scale (°F/°C)
- Water Resistant
- Large display
- Uses 2 "AAA" batteries (included)
- (Minimum 1000 full duty cycles over 24 month period)
- Temp range 95°F to 107.6°F +/- .2°F (35°C to 42°C +/- .1°C)
- Compact design measures just 4.7" x 1.17" x 1" (12 cm x 3.0 cm x 2.4 cm ) weighs 1.3oz (37g) including batteries



### ORDERING INFORMATION

#427	Adtemp™ Temple Touch Thermometer, Box
#427Q	Adtemp™ Temple Touch Thermometer, Display Pkg.

1-800-ADC-2670  
Visit us at [www.adctoday.com](http://www.adctoday.com)

