

Adtemp™ 424 Tympanic Thermometer

me

de

ADC

424

A

Large Backlit LCD Display



Audible "fever" alarm for temperatures above 99.5°F/37.5°C.



Switchable Dual scale Measures in either °F or °C modes.

Compact Ergonomic Design & Storage Caddy




Convenient Probe Cover Dispenser Contains 20 Covers

ADTEMP™ 424 INFRARED EAR THERMOMETER

Easy to use tympanic thermometer measures core temperature instantly. Convenient, economical probe covers make it the ideal instrument for use on all types of patients.

ADTEMP™ 424 Tympanic thermometer features:

- 1 second measurement using scan/peak technology
- Large 6mm illuminated display
- Temp range from 89.6°F - 108°F (32°C - 42.2°C)
- F/C switchable display
- Fever alarm at 99.5°F (37.5°C)
- Last reading memory
- Audible beep when measurement complete
- Auto off conserves battery strength
- Compact, ergonomic design - measures 5½" (13.97cm) long, weighs 1.6oz. (45.3g) including battery
- Powered by 1.55 (CR20320) Lithium battery (included)
- Includes 20 probe covers in convenient dispenser
- Desktop storage caddy



ORDERING INFORMATION

#424	Adtemp™ Infrared Ear Thermometer
#425-40	Adtemp™ Tympanic Sheaths, 40/pkg



THERMLIT120

1-800-ADC-2670
Visit us at www.adctoday.com

ADC
American Diagnostic Corporation

Adtemp™ 427 Temple Touch Thermometer

TM

Large Backlit
LCD Display

Compact
Ergonomic Design
(Actual Size)




ADTEMP™ 427 TEMPLE TOUCH THERMOMETER

Provides an accurate 6 second non-invasive temperature reading at the forehead making it ideal for pediatric patients. Cost effective design makes it ideal for institutional use.

ADTEMP™ 427 Temple Touch thermometer features:

- Non-invasive forehead temperature measurement using patented conductive technology delivers clinically accurate readings in about 6 seconds
- Reduces the risk of cross contamination and avoids the need for costly probe covers
- Ideal for the hospital, clinic, or physician office
- Auto off function conserves battery life
- Audible tone when measurement is complete
- Last reading memory
- Dual Scale (°F/°C)
- Water Resistant
- Large backlit display
- Uses 2 "AAA" batteries (included)
(Minimum 1000 full duty cycles over 24 month period)
- Temp range 95°F to 107.6°F (35°C to 42°C)
- Compact design measures just 4.7" x 1.17" x 1" (12 cm x 3.0 cm x 2.4 cm) weighs 1.3oz (37g) including batteries



ORDERING INFORMATION

#427	Adtemp™ Temple Touch Thermometer, Box
#427Q	Adtemp™ Temple Touch Thermometer, Display Pkg.

1-800-ADC-2670
Visit us at www.adctoday.com

