

Adtemp™ 429 Non-Contact Thermometer

Easy to use IR thermometer measures temperature instantly at the forehead!
Convenient, non-contact measurement makes the Adtemp™ 429 the ideal instrument for use on all patients.

ADTEMP™ 429 IR thermometer features:



- Reads in 1 second
- Large 6mm illuminated display
- Measuring range (Human)
50.0°F - 122.0°F +/- .5°F (10°C - 50°C +/- .3°C)
- Measuring range (Object)
32.0°F - 212.0°F +/- 1°F (0°C - 100°C +/- 1°C)
- °F/°C switchable display
- Switchable Human/Object modes
- 30 reading memory
- Night mode function mutes audibles for patient care
- Audible beep when measurement complete
- Auto off conserves battery strength
- Compact, ergonomic design - measures 5" (12.7cm) long,
- Weighs 3oz. (87g) including battery
- Powered by (2X) 1.55V AAA alkaline batteries (included)
- Optional wall bracket available



Measurement Sensor Lens




The Adtemp™ 429 offers 2 scan modes:



ORDERING INFORMATION

#429	Adtemp™ Non-Contact Thermometer
------	---------------------------------



6 34782 08381 6

429LIT123

1-800-ADC-2670
Visit us at www.adctoday.com

ADC
American Diagnostic Corporation