

# Satin™ Disposable Laryngoscopes

The first TRULY Cost Effective disposable laryngoscope blade AND handle system on the market.

**Fiberoptic Style Disposable Blades** in Macintosh and Miller profiles available in the most popular sizes. Compatible with any brand of disposable or reusable handle that meets ASTM F965 or ISO 7376-1. Polycarbonate design provides the rigidity needed with greater patient safety than all-metal or hybrid blades.

**LED Illuminated Disposable Handle** powered by 3 "AAA" batteries (included), provides brighter, whiter, cooler light than halogen. Compatible with any brand of disposable or reusable blade that meets ASTM F965 or ISO 7376-1.

Each blade and handle packaged in a hermetically sealed pouch ready for ETO sterilization. Packed 20 instruments per carton, 120 (6 cartons) per case.

**MSRP for any size blade or handle is JUST \$4.99 each (\$9.99 for 1 blade and 1 handle), making it the most economical disposable laryngoscope system on the market.**




1-800-ADC-2670  
Visit us at [www.adctoday.com](http://www.adctoday.com)

  
American Diagnostic Corporation

Satin™

# Satin™ Disposable Laryngoscopes



## MACINTOSH DISPOSABLE FIBEROPTIC

- Polycarbonate construction
- Fiberoptic light transmission provides cool, obstruction free viewing
- Featuring the 3 most popular curved blade sizes (2-4)
- Individually hermetically sealed, 20/carton



| ORDERING INFORMATION |       |                                |       |       |
|----------------------|-------|--------------------------------|-------|-------|
| Item#                | Blade | Description                    | Dim A | Dim B |
| 4072D-20             | 2     | Child                          | 114.3 | 89    |
| 4073D-20             | 3     | Med. Adult                     | 127   | 105   |
| 4074D-20             | 4     | Lrg. Adult                     | 155.5 | 133   |
| 4079D                | 2-4   | Kit 2 of each size with Handle |       |       |


(Dimensions in millimeters)

## MILLER DISPOSABLE FIBEROPTIC

- Polycarbonate construction
- Fiberoptic light transmission provides cool, obstruction free viewing
- Featuring the 4 most popular straight blade sizes (0-3)
- Individually hermetically sealed, 20/carton

| ORDERING INFORMATION |       |                                |       |       |
|----------------------|-------|--------------------------------|-------|-------|
| Item#                | Blade | Description                    | Dim A | Dim B |
| 4080D-20             | 0     | Preemie                        | 79.3  | 54    |
| 4081D-20             | 1     | Infant                         | 104.7 | 80    |
| 4082D-20             | 2     | Child                          | 152.4 | 133   |
| 4083D-20             | 3     | Adult                          | 196.8 | 171   |
| 4089D                | 0-3   | Kit 2 of each size with Handle |       |       |

(Dimensions in millimeters)



## LED DISPOSABLE BATTERY HANDLE

- LED illuminated for bright, cool light
- Ribbed surface for positive grip
- Sized like a traditional "C" cell handle
- Powered by 3 "AAA" batteries (included)
- May be ETO sterilized
- Each in hermetically sealed pouch, 20/carton
- The lowest cost handle on the market



| ORDERING INFORMATION |                               |
|----------------------|-------------------------------|
| Item#                | Description                   |
| 4065D-20             | LED Disposable Battery Handle |

1-800-ADC-2670  
Visit us at [www.adctoday.com](http://www.adctoday.com)



DISPOSLARYNGLIT