

# ADCUFF S.P.U.<sup>TM</sup>

Blood Pressure Cuffs



Single Patient Use  
Antibacterial  
Bladderless  
Latex Free



Adult and Pediatric ADCUFF S.P.U.™ in six clinical sizes and a variety of configurations to suit virtually any need, or budget.

#### ADCUFF S.P.U.™ features:

- Index and range marking to assure use of correct size
- Artery mark simplifies correct positioning
- Color coded sizes facilitate cuff selection
- Latex free construction reducing risk of allergic reaction
- Ultra-soft lining for patient comfort
- Individually leak tested to assure dependability
- Durable enough for repeated use, economical enough for single patient applications
- Vinyl material treated with antibacterial agent
- Six clinical sizes to choose from all in accordance with American Heart Association (AHA) recommendations
- One and two tube configuration for versatility
- A wide selection of connectors and adaptors for use with most blood pressure instruments
- Available in complete inflation system configurations
- Individually poly-bagged, 10/box
- Assembled in the U.S.A.

#### ADCUFF S.P.U.™ Ordering Chart

Size	Range	Color	Two Tube			One Tube			Inflation System
			Slip Connector (8972)	Screw Connector (8971)	Female Connector (8973)	Slip Connector (8972)	Female Connector (8973)	Bayonet Connector (8974)	Slip Connector (891F)
Infant	8-13cm	light blue/white	#8400I2-10	#8400I1-10	#8400I3-10	#8500I3-10	#8500I1-10	#8500I2-10	#8600I-10
Child	12-19cm	lavender/white	#8400C2-10	#8400C1-10	#8400C3-10	#8500C3-10	#8500C1-10	#8500C2-10	#8600C-10
Small Adult	17-25cm	green/white	#8400SA2-10	#8400SA1-10	#8400SA3-10	#8500SA3-10	#8500SA1-10	#8500SA2-10	#8600SA-10
Adult	23-33cm	navy/white	#8400A2-10	#8400A1-10	#8400A3-10	#8500A3-10	#8500A1-10	#8500A2-10	#8600A-10
Adult, Long	23-33cm	orange/white	#8400LA2-10	#8400LA1-10	#8400LA3-10	#8500LA3-10	#8500LA1-10	#8500LA2-10	#8600LA-10
Large Adult	31-40cm	pink/white	#8400X2-10	#8400X1-10	#8400X3-10	#8500X3-10	#8500X1-10	#8500X2-10	#8600X-10

Refer to page 4 for connector diagrams





Neonatal ADCUFF S.P.U.™ in five clinical sizes and a variety of configurations to suit virtually any need, or budget.

**Neonatal ADCUFF S.P.U.™ features:**

- Index and range marking to assure use of correct size
- Artery mark simplifies correct positioning
- Color coded sizes facilitate cuff selection
- Latex free construction reducing risk of allergic reaction
- Ultra-soft lining for patient comfort
- Individually leak tested to assure dependability
- Durable enough for repeated use, economical enough for single patient applications
- Vinyl material treated with antibacterial agent
- Five neonatal sizes to choose from all in accordance with American Heart Association (AHA) recommendations
- One and two tube configuration for versatility
- Equipped with male luer slip connectors - additional connectors available
- Individually polybagged, 10/box
- Assembled in the U.S.A.

**ADCUFF S.P.U.™ Ordering Chart**

Size	Range	Color	Two Tube	One Tube
			Slip Connector (8975)	Slip Connector (8975)
#1	3-6cm	light blue/white	#8401N-10	#8501N-10
#2	4-8cm	lavender/white	#8402N-10	#8502N-10
#3	6-11cm	green/white	#8403N-10	#8503N-10
#4	7-13cm	navy/white	#8404N-10	#8504N-10
#5	8-15cm	pink/white	#8405N-10	#8505N-10

Refer to page 4 for connector diagrams

## Accessories



### Latex Free Single Patient Use Bulb and Valve Assembly

#### The S.P.U. Bulb and valve features:

- Latex free PVC construction
- All metal end valve for dependability
- PVC deflation control valve
- 10/box
- Assembled in the U.S.A. from foreign components

#875PLF-10	Bulb and valve, S.P.U.
------------	------------------------



### ADSCOPE™ 665 Single Patient Use Diaphragm Only Stethoscope










#### The 665Y features:

- Diaphragm only chestpiece design
- Remarkable acoustic clarity in a lightweight, single patient use design
- 22" tubing
- Latex free construction
- Individually polybagged
- 10/box

#665Y-10	ADSCOPE™ 665Y, Stethoscope, S.P.U.
----------	------------------------------------



## Connectors and Adaptors

	Screw adaptor with barbed end for use with Dynamap <b>#8971-10</b>
	Male luer slip connector with barbed end (plastic) <b>#8972-10</b>
	Male luer slip connector with barbed end for use with ADC & Baum (metal) <b>#891M-10</b>
	Female locking luer connector with barbed end (plastic) <b>#8973-10</b>
	Female luer slip connector with barbed end for use with ADC & Baum (metal) <b>#891F-10</b>
	Bayonet connector for use with Hewlett Packard and Siemens monitors <b>#8974-10</b>
	Neonatal cuff connector, male luer slip fitting for 3/32" ID tubing <b>#8975-10</b>
	Double Male connector for 5/32" ID tubing (plastic) <b>#920-011-10</b>
	Y Adaptor, no fittings for use with two tube cuffs <b>#611-20-10</b>

55 Commerce Drive  
Hauppauge, NY 11788  
email: [info@adctoday.com](mailto:info@adctoday.com)  
<http://www.adctoday.com>  
NEW YORK • LONDON • TAIPEI • TOKYO