

Non-Contact Thermometer Cross Reference



| Brand/Model | ADC 429 | Exergen TAT-2000C | Graham Field 2220 | LinkTemp LMP001 | Mabis 18-535-000 | MDF NT13 | ThermoMedics Pro-TF300 | Veridian 09-349 | Welch Allyn 105801 |
|-------------------------|---------------------------------|------------------------------------|----------------------------------|----------------------|---------------------------------|--------------------------------|---------------------------|---------------------|-----------------------|
| Readout Time | 1 Second | Approx. 0.04 Seconds | "a few seconds" | 1 Second | 1 Second | 1 Second | 1-2 Seconds | 1 Second | 1 Second |
| Switchable F/C | Dual Scale | Dual Scale | Dual Scale | Dual Scale | Dual Scale | Dual Scale | Dual Scale | Dual Scale | Dual Scale |
| Modes | Body/Object | Temporal | Body/Object | Body/Object | Body/Object | Body/Object | Body/Object | Body/Object | Body |
| Proximity Sensor | Yes | No | No | No | Yes | No | No | Yes | Yes |
| Memory Recall | 30 Readings | 8 Readings | 10 Readings | 32 Readings | 10 Readings | 10 Readings | 30 Readings | 10 Readings | Last Reading |
| Fever Alarm | Yes | No | Yes | Yes | Yes | Yes | No | Yes | No |
| Display | Backlit | Backlit | N/A | Backlit | Backlit | Backlit | Backlit | Backlit | Backlit |
| Auto-Off | Yes | No | Yes | Yes | Yes | No | No | Yes | Yes |
| Battery | 2 AAA | 9V | 2 AAA | 2 AA | 2 AAA | 2 AAA | 2 AA | 2 AAA | 2 AAA |
| Battery Life | 2,000 Uses | 7,500 Uses | N/A | N/A | N/A | N/A | >15,000 Uses | N/A | 3,300 Uses |
| Warranty | 2 Years | 1 Year | 1 Year | 18 Months | 5 Years | 2 Years | 2 Years | 5 Years | 2 Years |
| Special Features | Night Mode/ Mutable Tones | Beeper/LED Can Be Turned Off | Retail Packaging/ Handheld | Lanyard | Retail Packaging Readings | Memory Can Average Stamp | Memory With Date/Time | Retail Packaging | Pulse Timer |
| Dimensions | 2.06" x 1.75" x 5" | 7" x 1.75" x 1.25" | 5.3" x 3.39" x 1.25" | 6.81" x 2.05" x 1.7" | 5.25" x 3.38" x 1.25" | N/A | 5.9" x 1.91" x 2.17" | 8.6" x 5.4" x 3" | 6.3" x 1.7" x 2.2" |
| Weight | 3 oz | 4.5 oz | 4.3 oz | 3.48 oz | 4.27 oz | N/A | N/A | 8 oz | 2.8 oz |
| MSRP | \$79.99 | \$99.99 | \$79.99 | \$170 | \$89.99 | \$79.99 | >\$300 | \$85.55 | >\$200 |