

Connect with Your Heart Health Goals

ADC's newest upper-arm digital home blood pressure monitor combines best-in-class features with Bluetooth connectivity!



Adult Wide-Range D-bar Cuff
22 to 42 cm. (8.7 to 16.5 in.)
Small Adult and Large Adult cuffs sold separately.

Displays
Systolic Pressure,
Diastolic Pressure and Pulse

The first blood pressure tracking app good enough to wear the ADC name.



ADC@Home Connectivity

Use Bluetooth to connect to our free app, where you can capture, graph, and store blood pressure readings and share with your healthcare providers.



adctoday.com/ADCatHome

Engineered for Comfort

Automatic inflation and deflation control with measurement performed on inflation for more comfortable readings.

Two-User Memory

Advanced memory with date and time stamp recalls 99 readings for two separate users -- perfect for couples!

Color-Coded Display

Using a built-in, color-coded traffic light indicator, this advanced monitor visually compares readings to standards set by the American Heart Association®.

Integrated, One-Button Average Mode

A single measurement is rarely your real result. The American Heart Association recommends two or more readings with the average being your real result. Advantage Averaging Mode takes three consecutive measurements, and calculates the average.

Irregular Heartbeat Detection

Detects the appearance of an irregular heartbeat during measurement and gives a warning signal with your reading. An irregular heartbeat may cause or indicate cardiovascular disease, so early detection is important.

For additional information, contact your ADC sales representative or our Customer Service team at 1-800-ADC-2670.

