

REPLACEMENT CUFF FOR ADC ADVANTAGE HOME BLOOD PRESSURE MONITORS

This cuff is designed for and ready to use with these ADC home blood pressure monitors:

- *Advantage 6112*
- *Advantage Plus 6121*
- *Advantage Connect 6024N*

To connect the cuff, simply insert the connector into the hole in the left side of the monitor. Make sure the cuff connector is fully inserted before taking a reading.

This cuff is compatible with the legacy ADC home blood pressure models listed below, but you must change the connector before using:

- *Advantage Semi-Automatic 6012N*
- *Advantage 6021N*
- *Advantage Plus 6022N*

Please contact ADC Customer Service at 1-800-ADC-2670 or email info@adctoday.com for more information.