



ADC Advantage Home Blood Pressure Line for 2026



6115
Advantage Connect®



6112
Advantage™



6121
Advantage™ Plus



6024N
Advantage Connect®

Replaces	Models 6015 and 6016	Model 6012	Models 6021 and 6022	
Type	Auto Wrist	Auto Upper Arm	Auto Upper Arm	Auto Upper Arm
Measures Systolic/Diastolic/Pulse	Yes	Yes	Yes	Yes
Audible Readings	Yes	No	No	No
Display	LCD with Backlight	LCD	LCD with Backlight	LCD
Display Size	1.7" × 1.6" (43 mm × 40 mm)	2.6" × 2.5" (66.2 mm × 62.5 mm)	2.6" × 2.5" (66.2 mm × 62.5 mm)	2.75" × 2.5" (70 mm × 63.5 mm)
Systolic/Diastolic Range	60-230 mmHg / 40-130 mmHg	60-255 mmHg / 40-200 mmHg	60-255 mmHg / 40-200 mmHg	60-255 mmHg / 40-200 mmHg
Pulse Range	40-199 BPM	40-199 BPM	40-199 BPM	40-199 BPM
Pressure Accuracy	±3 mmHg	±3 mmHg	±3 mmHg	±3 mmHg
Pulse Accuracy	±5%	±5%	±5%	±5%
Power Source	2 AAA Batteries (included)	4 AAA Batteries (included)	4 AAA Batteries (included)	4 AAA Batteries (included)
USB-C Power Port	No	No	Yes	No
Memory	Last 199 Readings x2	Last 99 Readings x2	Last 99 Readings x2	Last 99 Readings x2
Irregular Heartbeat	Yes	Yes	Yes	Yes
MAM (Average Mode)	Yes	No	Yes	Yes
Traffic Light Indicator	Yes	No	Yes	Yes
Bluetooth	Yes	No	No	Yes
ADC@Home App	Yes	No	No	Yes
Auto-off / Sleep Mode	Yes	Yes	Yes	Yes
Dimensions	2.4" × 3.4" × .75" (6 cm × 8.6 cm × 2 cm)	5.25" × 3.65" × 2" (13.34 cm × 9.3 cm × 5 cm)	5.25" × 3.65" × 2" (13.34 cm × 9.3 cm × 5 cm)	5.35" × 3.5" × 2.4" (13.6 cm × 9 cm × 6 cm)
Weight	.23 lbs (106 g), excluding batteries	.83 lbs (378 g), including batteries and cuff	.86 lbs (390 g), including batteries and cuff	.60 lbs (270.5 g), including batteries and cuff
Latex Content	Not made with natural rubber latex	Not made with natural rubber latex	Not made with natural rubber latex	Not made with natural rubber latex
Cuff Type	Wrist	Comfort Contour™ with Self-Adjusting D-Bar	Comfort Contour™ with Self-Adjusting D-Bar	D-Bar, Self-Adjusting
Standard Cuff Size	5.5" to 8.5" (13.5 cm to 21.5 cm)	Adult, Extended Range: 8.7" to 16.5" (22 cm to 42 cm)	Three Standard Options: 6121 – Adult, Extended Range: 8.7" to 16.5" (22 cm to 42 cm) 6121X – Large Adult: 12.5" to 20.5" (32 cm to 52 cm) 6121SA – Small Adult: 6.7" to 8.7" (17 cm to 22 cm)	Adult, Extended Range: 8.7" to 16.5" (22 cm to 42 cm)
Additional Cuff Sizes Available	No	Small Adult (#850-6121SA) Large Adult (#850-6121X) Sold Separately	See Above	Small Adult (#850-6121SA) Large Adult (#850-6121X) Sold Separately
Storage Case	No	No	No	No
Warranty	3 Year	5 Year	5 Year	5 Year
MSRP	\$54.99	\$49.99	\$59.99	\$93.04

Contact your ADC sales representative at **1-800-ADC-2670** or visit **www.adctoday.com/ADCatHome**