

One of the
**Most Advanced,
One-piece Blood Pressure Cuffs**
On the Market



ADCUFF™ +

Our all-new line of one-piece, reusable blood pressure cuffs is designed for measurement accuracy, patient comfort, and practitioner ease.

- Simple to use and easy to clean
- Bladderless construction from TPU-coated, 70 denier nylon
- Surface-mounted ports for a uniform fit around patient limb
- Longer limb ranges for improved patient coverage
- Hang tab (on most sizes) fits conventional pocket gauges
- Hook-and-loop adhesive rated for tens of thousands of cycles
- Not made with natural rubber latex or phthalates
- Three-year warranty

Nylon construction can be easily cleaned with most major brands of alcohol- or bleach-based wipes. Minimize risks of cross contamination

The perfect solution when versatility and ease of cleaning is critical



www.adctoday.com/AdcuffPlus

ADC
AMERICAN DIAGNOSTIC CORPORATION

ADCUFF™+

Designed for measurement accuracy,
ease of use, and patient comfort


Durable Hook & Loop

Withstands tens of
thousands of
open/close cycles

Curved Corners
Improves patient
comfort

**Industry-Standard
Color Coding**
Simplifies stock control
and size identification

Gauge Hang Tab
For use with traditional
pocket aneroids



Cuff ID Panel
Provides vital cuff information
at a glance

Part Number Reference
Simplifies stock control
and reordering

Surface-Mounted Tube Ports
Improves durability and allows
a better patient fit

Four-Sided Welding
Lies flatter
for a better fit

Also Featuring ADC's Size Guide™ Marking System



Size icons helps ensure
proper cuff size and range
selection



**Properly positioned
artery mark** – you'd
be surprised how
many aren't



**Accurate, comprehensive
range markings**, helps
prevent miscuffing

Convenient and convertible

- Six clinical cuff sizes:
Infant, Child, Small Adult,
Adult, Large Adult,
and Thigh
- Single-tube and double-
tube designs
- Complete inflation
systems
- Adcuff+ versions of many
ADC sphyg models

Visit
www.adctoday.com/AdcuffPlus
to view all model options



700-13TBR1



705PMK-MCC1



750D-11AN1



728-619-11AN1

For additional information, contact your ADC sales representative
or our Customer Service team at 1-800-ADC-2670.

