

ADView[®] 2.6

Clinical-Grade Diagnostic
Station from ADC[®]

NEW!



Desktop

The Advview 2.6 features a modern, compact design that's self-standing for desktop or tabletop use.

Mobile

Optimize your workflow with our rugged, easy to assemble, five-legged mobile stand. It's height-adjustable with an oversized basket to keep accessories organized and features large wheels – two with locking casters.



Wall Mount

Choose the custom mounting plate with integral basket for easy wall mounting.

PRACTICAL INNOVATION FOR PATIENT CARE

See critical patient data at a glance. Big, bright, 5.6" high-resolution color LCD display readings easy to view in a variety of clinical environments.

Go wherever patient care happens. Lightweight, compact construction makes it easy to move from room to room or between departments.

Spend less time navigating. An intuitive interface helps clinicians access measurements quickly with minimal training.

Fits the way you work. Flexible configuration options allow the monitor to adapt to a wide range of care settings and workflows.

Measure blood pressure your way. Multiple measurement modes provide the flexibility to support routine spot checks and multi-reading averages.

Track changes as they happen. View current measurements alongside trends to help identify changes in patient condition over time.

Power through an entire shift. A long-lasting rechargeable Li-ion battery delivers up to 14 hours of continuous operation, reducing downtime and supporting mobile care.

One monitor for virtually every patient. Designed to support adult and pediatric populations, helping simplify equipment needs across the care continuum.

Contact your **ADC Sales Representative**
or our **Customer Service Team** at
1-800-ADC-2670

adctoday.com/adview26





BLOOD PRESSURE • TEMPERATURE • SPO2 • PULSE RATE



9006BPS
Blood Pressure & SpO2



9006BPSTO
Blood Pressure, SpO2 & Oral/Axillary Thermometry



9006BPSTT
Blood Pressure, SpO2 & IR Tympanic Thermometry

ORDERING INFORMATION

ADVIEW 2.6 DIAGNOSTIC STATION

BP, SpO2	9006BPS
BP, SpO2, Oral/Axillary Therm	9006BPSTO
BP, SpO2, Infrared Tympanic Therm	9006BPSTT
BP, SpO2, Printer	9006BPSP
BP, SpO2, Oral/Axillary Therm, Printer	9006BPSTOP
BP, SpO2, Infrared Tympanic Therm, Printer	9006BPSTTP

BP base unit includes: Rechargeable lithium-ion battery; U.S. power cord; 10-foot length of tubing; Adult and Large Adult Adcuff+ cuffs

Oral/Axillary Temperature module includes: Probe, well, and box of 25 probe covers

IR Temperature module includes: Thermometer, cradle, probe cover dispenser and box of 40 probe covers

SpO2 module includes: Adult finger sensor and 8-foot extension cable

MOUNTING PLATFORMS

Wall Mount with Basket	9006WB
Mobile Stand with Basket	9006M2

CUFF LIMB RANGES

Child 13-19.5 cm 9005-9CGR-1MB Sold Separately	Small Adult 19-27 cm 9005-10SARB-1MB Sold Separately	Adult 23-40 cm 9005-11AN-1MB Adult and Large Adult Cuffs Included	Lg. Adult 34-50 cm 9005-12XBD-1MB Sold Separately	Thigh 40-66 cm 9005-13TBR-1MB Sold Separately
--	--	---	---	---

TECHNICAL SPECIFICATIONS

Dimensions:	9.25" (L) × 6.85" (W) × 7.44" (H) (235 mm × 174 mm × 189 mm)
Weight:	Approx. 5.5lbs (2.49 kg) not including battery
BP Range:	Adult mode Dia: 30-215 mmHg / Sys: 60-270 mmHg
BP Accuracy:	Meets or exceeds ANSI/AAMI/ISO 81060-2:2013 for non invasive accuracy (±5 mmHg mean error with 8 mmHg SD)
Pulse Rate (SpO2)	
Range:	25 bpm - 300 bpm
Accuracy:	±2 bpm
Pulse Rate (NIBP)	
Range:	40 bpm - 240 bpm
Accuracy:	±3 bpm or 3.5%, whichever is greater
Temperature (Oral/Axillary)	
Range:	77°F ~ 113°F (25°C ~ 45°C)
Accuracy:	77°F ~ 113°F ± 0.2°F (25°C ~ 45°C ± 0.1°C)
Temperature (Infrared)	
Range:	93.2°F ~ 108°F (34°C ~ 42.2°C)
Accuracy:	95°F ~ 107.6°F ± 0.4°F (35.5°C ~ 42°C ± 0.2°C)
SpO2	
Range:	0%-100%
Accuracy:	±2% (70% to 100%)

WARRANTY

Main Unit, SpO2 modules, Temp module	3 Years
Blood Pressure Cuffs	3 Years
Wall Mount, Mobile Stand	2 Years
Accessories (SpO2 sensor, printer)	90 Days

Additional specifications online at
adctoday.com/products/adview26