

# advantage connect<sup>®</sup> 6024N

## Automatic Digital Home Blood Pressure Monitor

### FULL-FEATURED ACCURACY & CONNECTED CONVENIENCE

The **Advantage Connect<sup>®</sup> 6024N** delivers clinically precise blood pressure readings with smart, connected convenience. This upper-arm monitor features MAM Advantage™ triple-measurement averaging for enhanced accuracy, along with a color-coded hypertension indicator for at-a-glance insights. Bluetooth<sup>®</sup> connectivity syncs seamlessly with the free ADC@Home app, allowing users to record, track, and securely share readings over time. Comfortable, intuitive, and feature-rich, the 6024N is an ideal choice for users who want dependable insights and greater confidence in managing their heart health.

Updated for 2026



5-Year  
Warranty

Adult Extended-Range cuff fits arm circumferences from 8.7" to 16.5" (22 cm to 42 cm)  
Small Adult and Large Adult cuffs sold separately

#### ✓ Measures on Inflation

Automatic inflation and deflation control with measurement performed on inflation for more comfortable readings.

#### ✓ MAM Advantage™ Averaging Mode

A technique recommended by the American Heart Association<sup>®</sup> for enhanced accuracy, our automated averaging mode takes three consecutive measurements then calculates and displays a single result.

#### ✓ Two-User Memory

Stores up to 99 readings per user with date and time stamps, making it easy to track trends over time.

#### ✓ Traffic Light Indicator

This built-in color-coded scale visually aligns readings with the American Heart Association<sup>®</sup> guidelines for quick interpretation.

#### ✓ Irregular Heartbeat Detection

Alerts you to potential irregularities so you can share meaningful information with your healthcare provider.

### Premium Comfort Contour™ Cuffs

At the heart of our upper-arm monitors is the thoughtfully engineered Comfort Contour™ cuff — crafted from soft, premium nylon and shaped to follow the natural curve of your arm. Unlike flat or straight cuffs, its ergonomic contoured design distributes pressure evenly along the arm, eliminating pinching and bunching while ensuring the bladder inflates consistently for a secure fit and more accurate readings.

### ADC@Home™ Connectivity

Connect via Bluetooth<sup>®</sup> to our free app, where you can capture, graph, and store blood pressure readings and share with your healthcare providers. Learn more at [www.adctoday.com/ADC@Home](http://www.adctoday.com/ADC@Home).



Contact your ADC sales representative at **1-800-ADC-2670** or visit [www.adctoday.com/products/6024N](http://www.adctoday.com/products/6024N)