

advantage™ 6112

Automatic Digital Home Blood Pressure Monitor

RELIABLE EVERYDAY MONITORING

The **Advantage™ 6112** is your reliable partner for proactive health management. Designed for simplicity, comfort, and confidence, it combines advanced measurement technology with intuitive operation. Whether you're just beginning to monitor your blood pressure or maintaining an ongoing routine, the 6112 offers the features you need at an affordable price point.



Extra-large LCD screen measures 2.6" x 2.5"

Systolic pressure

Diastolic pressure

Pulse

4 AAA batteries (included)



5-Year Warranty

Adult, Extended Range cuff fits arm circumferences from 8.7" to 16.5" (22 to 42 cm)
Small Adult and Large Adult cuffs sold separately

✓ Measures on Inflation

Automatic inflation and deflation control with measurement performed on inflation for more comfortable readings.

✓ Automatic Fit Check

Built-in fit detection ensures your cuff is positioned correctly before each reading, reducing user error.

✓ Two-User Memory

Stores up to 99 readings per user with date and time stamps, making it easy to track trends over time.

✓ Irregular Heartbeat Detection

Alerts you to potential irregularities so you can share meaningful information with your healthcare provider.

✓ Easy-to-Read Display

Large, high-contrast screen with oversized numerals simplifies viewing for all users.



Premium Comfort Contour™ Cuffs

At the heart of our upper-arm monitors is our thoughtfully engineered Comfort Contour™ cuff, crafted from soft, high-quality nylon and ergonomically contoured to fit the natural shape of your arm. Unlike flat, straight cuffs, a curved or conical cuff helps distribute pressure more evenly along the arm's contour, reducing pinching or bunching and ensuring the bladder inflates consistently over the arm for a better fit and more accurate readings.



Contact your ADC sales representative at **1-800-ADC-2670** or visit www.adctoday.com/products/6112