


# advantage connect<sup>®</sup> 6115

## Automatic Digital Wrist Home Blood Pressure Monitor

### LIGHTWEIGHT PORTABILITY & CONNECTED TRACKING

Compact enough to fit in the palm of your hand, the **Advantage Connect<sup>®</sup> 6115** wrist monitor brings together smart technology and clinical precision for effortless, on-the-go blood pressure monitoring. Lightweight, connected, and easy to use, it's the perfect choice for users who want reliable results wherever life takes them.



4 AAA batteries (included)

Attached wrist cuff fits arm circumferences from 5 1/3" to 8 1/2" (13.5 cm to 21.5 cm)

#### ✓ Measures on Inflation

Automatic inflation and deflation control with measurement performed on inflation for more comfortable readings.

#### ✓ Automatic Fit Check

Built-in fit detection ensures your cuff is positioned correctly before each reading, reducing user error.

#### ✓ MAM Advantage™ Averaging Mode

A technique recommended by the American Heart Association<sup>®</sup> for enhanced accuracy, our automated averaging mode takes three consecutive measurements then calculates and displays a single, result.

#### ✓ Two-User Memory

Stores up to 199 readings per user with date and time stamps, making it easy to track trends over time.

#### ✓ Guest Mode

A convenient option for single readings that don't need to be stored in memory.

#### ✓ Traffic Light Indicator

This built-in color-coded scale visually aligns readings with the American Heart Association<sup>®</sup> guidelines for quick interpretation.

#### ✓ Irregular Heartbeat Detection

Alerts you to potential irregularities so you can share meaningful information with your healthcare provider.

#### ✓ Easy-to-Read Display

Large, high-contrast screen with oversized numerals simplifies viewing for all users.

### ADC@Home™ Connectivity

Connect via Bluetooth<sup>®</sup> to our free app, where you can capture, graph, and store blood pressure readings and share with your healthcare providers. Learn more at [www.adctoday.com/ADC@Home](http://www.adctoday.com/ADC@Home).



Contact your ADC sales representative at **1-800-ADC-2670** or visit [www.adctoday.com/products/6115](http://www.adctoday.com/products/6115)