

# Comparison of ADC's PMV Otoscope to WA MacroView™ and our Standard Otoscope



FEATURES	ADC Standard Otoscope	ADC PMV Otoscope	WA MacroView™ Otoscope
Field of View	5.5 mm	11 mm	10.5 mm
Magnification	3X	5X	4X
Technology	Simple Lens	Prism	Compound Lens
Focusing Wheel	N/A	Yes	Yes
Interchangeable Bayonet Mount	Yes	Yes	Yes
Insufflator Connection	Barb	Barb	Port
Lighting Options	LED/Hal-Xenon	LED/Hal-Xenon	LED/Hal-Xenon
Instrument Warranty	2 Year	2 Year	2 Year
LED Warranty	Life	Life	10 Year