

# ***e-sphyg3™***

Digital Blood Pressure Monitor

- **Accurate**
- **Affordable**
- **Clinical Grade**



ADC's all-new e-sphyg 3 captures blood pressure measurements with the touch of a button. Designed for busy professionals who value both accuracy and affordability, the e-sphyg 3 is the feature-packed, clinical-grade solution you've been waiting for.

### Three clinically validated modes:



#### Single measurement mode

with pulse rate and mean arterial pressure (MAP) displayed in as little as 30 seconds.



#### Customizable, three-reading average mode

with individual measurement display.



#### Manual mode

for auscultation with a stethoscope.




Last-reading memory




Rechargeable high-capacity battery  
**STANDARD**



The e-sphyg<sup>3</sup> is perfect for:



Physician practices



Walk-in clinics



Ambulatory surgery centers

**Desktop design with built-in handle**  
Or choose one of these optional platforms  
to fit your practice workflow.

**Oscillometric measurements capture:**

**[SYS]** Systolic blood pressure  
(60 to 255 mmHg)

**[DIA]** Diastolic blood pressure  
(30 to 200 mmHg)

**[HR]** Pulse rate  
(40 to 200 BPM)

**[MAP]** Mean arterial pressure  
(35 to 240 mmHg)

**[IB]** Irregular  
Heartbeat

**For adult and pediatric patients**

**INCLUDES:** Three-cuff kit  
(Small Adult, Adult, and Large Adult) for accurate readings on virtually all patients.

Choose our traditional Adcuff design or our new Adcuff+ one-piece cuffs.

**Wall Mount**

Custom bracket for easy mounting to the wall or on an ADC Adstation system.

*Includes cuff storage basket.*



**Mobile Stand**

Our rugged, easy-to-assemble mobile stand features five rolling casters (two locking) for room-to-room portability.

*Includes cuff storage basket.*



**Desk Caddy**

A modern, compact design that helps keep cuffs and cords organized.



Dentist offices



Birth centers



Long-term care facilities



Digital Blood Pressure Monitor

The clinical-grade solution you've been waiting for.

#### Technical Comparisons

| FEATURE                          | e-sphyg3 | Vital Signs Monitors | Home BP Devices | Manual BP Devices |
|----------------------------------|----------|----------------------|-----------------|-------------------|
| Automated Measurement            | Yes      | Yes                  | Yes             | No                |
| Clinical Grade NIBP              | Yes      | Yes                  | No              | N/A               |
| 3 Reading Average Mode           | Yes      | Some                 | Some            | N/A               |
| Automatic and Manual Modes       | Yes      | Some                 | No              | No                |
| High-Capacity Battery - STANDARD | Yes      | No                   | No              | N/A               |
| Includes 3 cuffs                 | Yes      | No                   | No              | No                |
| Multiple Platforms               | Yes      | Yes                  | No              | Yes               |
| Irregular Heartbeat              | Yes      | Some                 | Some            | No                |



#### Ordering Information

**9003K-MCC** e-sphyg 3 with two-piece cuffs

**9003K-MCC1** e-sphyg 3 with one-piece cuffs

**Both configurations include:** e-sphyg 3 monitor; built-in, rechargeable, high-capacity battery; 6-foot power cord with four adapters; three Adcuff blood pressure cuffs (small adult, adult, large adult); and 4.5-foot tubing with connector.

#### Optional Platforms

**9003D** e-sphyg 3 desk caddy

**9003W** e-sphyg 3 wall mount bracket (compatible with ADC Adstation systems) and cuff basket

**9003M** e-sphyg 3 mobile stand and cuff basket

#### Replacement Parts

**9003ADAPT** AC adapter set for esphyg 3

**9003BAT** High-capacity battery for esphyg 3

**9003CK** Adcuff kit (small adult, adult, large adult) with e-sphyg 3 connectors

**9003CT** Extension tubing/connectors for esphyg 3



55 Commerce Drive  
Hauppauge, NY 11788

Toll Free: 1-800-ADC-2670  
(1-800 232-2670)  
Email: [info@adctoday.com](mailto:info@adctoday.com)  
Web: [www.adctoday.com](http://www.adctoday.com)